sebia

**CAPICLEAN CAPILLARYS 3
CAPILLARYS 3 CAPICLEAN**

Ref. 2060

IVD

CE

2019/12

Intended use

For the sample probe cleaning in automated instrument for capillary electrophoresis CAPILLARYS 3, SEBIA, during the CAPICLEAN cleaning sequence.

For *In Vitro* diagnostic use.

Composition / Reagent

ITEMS	PN 2060
CAPICLEAN (ready to use)	1 vial, 25 mL
Pierceable cap for the CAPICLEAN vial (cf "PHOTOS")	1 cap

Concentrated CAPICLEAN solution contains : proteolytic enzymes, surfactants and additives nonhazardous at concentrations used, necessary for optimum performances.

WARNING : The CAPICLEAN solution may cause irritation or burns to skin, eyes and mucous membranes. See the safety data sheet.

Storage / Stability

Store CAPICLEAN refrigerated (2 – 8 °C). It is stable until the expiration date indicated on the vial label. DO NOT FREEZE. Precipitate or combined particles in suspension (flocules) may be observed in the CAPICLEAN vial without any adverse effects on its utilization. Do not dissolve this precipitate or these particles. The CAPICLEAN vial placed in the secondary compartment of the CAPILLARYS 3 instrument is stable for a maximum of 3 months.

NOTE : During transportation, CAPICLEAN can be kept without refrigeration (15 to 30 °C) for 15 days without any adverse effects on its efficiency.


Procedure

See the instruction manual of CAPILLARYS 3, SEBIA.

1. Replace the cap of the CAPICLEAN vial by the pierceable cap. Check that the vial is perfectly closed (once in place, the pierceable cap may stay all the time on the CAPICLEAN vial if the vial is in a vertical position, do not turn the vial upside down to avoid leakages).


CAPICLEAN vial (front side)

CAPICLEAN vial (back side)

- Place the CAPICLEAN vial in the temperature controlled area of the CAPILLARYS 3 instrument secondary compartment in position S1, S2 or S3, place the side of the vial with the RFID label towards the left.


Position of the CAPICLEAN vial in the secondary compartment of CAPILLARYS 3


- Close the lid of the secondary compartment.
 - To start the CAPICLEAN cleaning cycle :
- Press "Maintenance & Cleaning" from the main menu :

- In the "Maintenance & Cleaning" menu, press "Start CAPICLEAN" to start the sample probe cleaning with the CAPICLEAN reagent.

- Confirm the startup of the cycle by pressing (= Yes), and if not, press (= No).

NOTES :

- the CAPICLEAN cycle will be performed when the instrument is waiting for rack and without any migration,
- the CAPICLEAN cycle lasts 10 minutes.

WARNING : In the current version of the instrument, the absence of progress bar or message indicating that the CAPICLEAN cycle is in progress may suggest that the instrument is inactive. Wait for the end of the CAPICLEAN cycle before any other processing.

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WARNING : The pierceable cap is composed of 2 parts : a screw cap and a septum. When the septum is separated from the screw cap, replace it in the cap as follows and check that it is correctly inserted before placing the pierceable cap on the reagent vial.

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sebia Swiss GmbH
sebia UK Ltd