

sebia

CLEAN PROTECT

Ref. 2059

IVD

CE

2019/12

For *In Vitro* Diagnostic Use.

Composition

The vial of CLEAN PROTECT concentrated solution contains a preservative.

WARNING : *The CLEAN PROTECT solution presents a characteristic odor. See the safety data sheets.*

Use

WARNING : *When dispensing the CLEAN PROTECT solution, wear protective work clothing, gloves and goggles.*

With CAPILLARYS and MINICAP systems

For the storage of distilled or deionized water intended for rinsing capillaries in automated systems CAPILLARYS and MINICAP, SEBIA, for capillary electrophoresis.

Add 400 μL of CLEAN PROTECT in one litre of distilled or deionized water (or 100 μL of CLEAN PROTECT in 250 mL of distilled or deionized water, 300 μL in 750 mL or 500 μL in 1.2 L). Mix well.

NOTE : It is recommended to filter distilled or deionized water with 0.45 μm filter before use.

WARNING : *Do not use marketed deionized water, such as water for ironing for example (risk of important capillaries damage). Use only water with ultrapure quality, such as injection grade water.*

Before filling the rinse container, it is recommended to wash it with plenty of distilled or deionized water.

With HYDRAGEL techniques

For the extended storage of the diluted destaining solution (for more than one week) to prevent microbial proliferation.

Add 250 μL of CLEAN PROTECT in 5 litres of diluted destaining solution (or 50 μL of CLEAN PROTECT in one litre of diluted destaining solution). Mix well.

For the sponge of HYDRASYS, HYDRASYS 2, CAPILLARYS, HYDRAPLUS and ASSIST humidity chamber

For the storage of the sponge of humidity chamber after a washing with plenty of water.

Add 90 μL of CLEAN PROTECT in 250 mL of distilled or deionized water. Mix well.

Apply the volume of diluted CLEAN PROTECT solution on the sponge of the humidity chamber (in place of the sodium azide solution) indicated in the instructions for use of each humidity chamber.

Storage, stability and signs of deterioration

Store CLEAN PROTECT at room temperature (15 – 30 °C) or refrigerated (2 – 8 °C). It is stable until the expiration date indicated on the vial label. DO NOT FREEZE.

CLEAN PROTECT must be free of precipitate. Discard CLEAN PROTECT if it changes its appearance.

For CAPILLARYS and MINICAP, distilled or deionized water with CLEAN PROTECT is stable for 1 month maximum at room temperature, in automated system or in closed container. Discard water if it changes its appearance, e.g., becomes cloudy due to microbial contamination.

On the sponge of humidity chambers, the diluted CLEAN PROTECT solution is stable for 1 month maximum at 2 – 8 °C. Change the solution after having rinsed the sponge with plenty of water according to the procedure described in the humidity chamber package insert.



Parc Technologique Léonard de Vinci
CP 8010 Lisses - 91008 EVRY Cedex - France
Tél. : 33 (0)1 69 89 80 80 - e-mail : sebia@sebia.com

sebia Benelux SCS / Comm. V

Jan Olieslagerstraat, 41
1800 Vilvoorde
Belgique / België
Tél. : 32 (0)2 702 64 64
Fax : 32 (0)2 702 64 60
e-mail : sebia.benelux@sebia.be

sebia Brasil.

Rua Barão do Triunfo, 73, Cj 74
CEP 04602-000
São Paulo
Brasil
Tel. : 55 11 3849 0148
Fax : 55 11 3841 9816
e-mail : sebia@sebia.com.br

sebia GmbH

Münsterfeldallee, 6
36041 Fulda
Deutschland
Tel. : 49 (0)661 3 30 81
Fax : 49 (0)661 3 18 81
e-mail : sebia@sebia.de

sebia Hispania s.a.

C/Sicilia, nº 394
08025 Barcelona
España
Tel. : 34 93 208 15 52
Fax : 34 93 458 55 86
e-mail : sebia@sebia.es

sebia Inc.

400-1705 Corporate Drive
Norcross, GA 30093
U.S.A.
Tel. : 1 770 446 - 3707
Fax : 1 770 446 - 8511
e-mail : info@sebia-usa.com

sebia Italia S.r.l.

Via Antonio Meucci, 15/A
50012 Bagno a Ripoli (FI)
Italia
Tel. : 39 055 24851
Fax : 39 055 0982083
e-mail : info@sebia.it

sebia Swiss GmbH

Verenastrasse, 4b
CH-8832 Wollerau
Switzerland
Tel. : 41 44 787 88 10
Fax : 41 44 787 88 19
e-mail : contact.ch@sebia.com

sebia UK Ltd

River Court, The Meadows Business Park
Station Approach, Blackwater
Camberley, Surrey, GU17 9AB
United Kingdom
Tel. : 44 (0)1276 600636
Fax : 44 (0)1276 38827
e-mail : sales@sebia.co.uk

sebia

Shanghai Representative Office
Cross Tower, Room 2306-07
318 Fuzhou Road
Shanghai 200001
China
Tel. : 00 86 (21) 6350 1655
Fax : 00 86 (21) 6361 2011
e-mail : sebia@sebia.cn