RIDASCREEN[®] FAST Ei / Egg Protein

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RIDASCREEN[®]FAST Ei / Egg Protein

Brief information

RIDASCREEN[®]FAST Ei / Egg Protein (Art. No. R6402) is a sandwich enzyme immunoassay for the quantitative analysis of whole egg powder in food.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction and centrifugation
Time requirement:	sample preparation (for 10 samples)..... approx. 20 min test implementation (incubation time)..... 30 min
Standard material:	NIST Reference Material 8445 (whole egg powder)
Limit of detection:	0.10 mg/kg (ppm) whole egg powder (corresponds to 0.03 mg/kg (ppm) egg white protein)
Limit of quantification:	0.5 mg/kg (ppm) whole egg powder (corresponds to 0.13 mg/kg (ppm) egg white protein)
Specificity:	The specific antibodies detect from the egg white proteins ovalbumin and ovomucoid. No cross reactivity was observed to raw or cooked meat of chicken, turkey hen, pork and beef.

1. Intended use

RIDASCREEN®FAST Ei / Egg Protein is a sandwich enzyme immunoassay for the quantitative analysis of whole egg powder in food like noodles, salad dressings, sausages, wines, baking-mixtures for cakes or bread and ice cream.

2. General

Egg can be present as an ingredient or as a contamination in raw and cooked products. According to the **EU Directive 2003/89/EG** from 10th November 2003, egg must be declared on food labels as it can induce allergic reactions.

Egg white contains 9 - 11 % of protein. Four main allergenic proteins come to 80 % of the egg white proteins. These main allergens are ovomucoid (11 %), ovalbumin (54 %), ovotransferrin (12 %) and lysozyme (3.5 %). The allergenic potential of the proteins present in egg yolk is only moderate.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to egg white proteins. By adding standards and samples to the wells, any egg white protein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of egg white protein takes place by adding Substrate/Chromogen solution. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the egg white proteins concentration of the sample. The result is expressed in mg/kg whole egg powder.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with antibodies against egg white proteins
- 5 x Standard solutions *) (1.3 ml each)
0 mg/kg (zero standard), 0.5 mg/kg, 1.5 mg/kg, 4.5 mg/kg, 13.5 mg/kg whole
egg powder in aqueous solution, ready to use

- 1 x Conjugate (0.7 ml)red cap
peroxidase conjugated antibody, 11fold concentrate
- 1 x Substrate/Chromogen solution (10 ml)..... brown cap
stained red
- 1 x Stop solution (14 ml).....yellow cap
contains 1 N sulfuric acid
- 1 x Allergen Extraction buffer (125 ml)green cap
as a **10fold concentrate**
- 1 x Washing buffer (100 ml) brown cap
as a 10fold concentrate

*) The **dilution factor of 20** for the sample has already been taken into account. Therefore, the whole egg powder concentration of samples can directly be read from the standard curve.

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal glass vials
- shaker
- water bath
- laboratory mincer/grinder, mortar, ultra-turrax or homogenizer
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl -micropipettes

5.2. Reagents:

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The washing buffer contains thimerosal. Avoid skin contact.

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate-/chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of egg and to avoid contamination.

The **Allergen Extraction buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99°F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water, or alternatively use the whole content of the bottle and fill up to 1250 ml). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. twelve weeks.

- grind 5 g of the sample carefully and mix thoroughly
- weigh 1 g of the sample and add 20 ml diluted Allergen Extraction buffer (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- or add 19 ml of the diluted Allergen Extraction buffer to 1 ml liquid sample (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- mix intensively and incubate for 10 min at 60 °C (140 °F), afterwards cool down
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl per well in the assay

Remark:

For cloves, mustard or celery containing samples, add 1 g casein to 1 g or 1ml sample to avoid disturbing effects.

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. The extracts can be stored at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **antibody enzyme conjugate** (bottle with red cap) is provided as an 11fold concentrate. Since the diluted enzyme conjugate solution has a limited stability, only the amount which actually is needed should be reconstituted. Before diluting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **washing buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted washing buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard solution or prepared sample to the wells, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete

removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

4. Add 100 µl of the diluted enzyme conjugate solution (see 10.1.) to each well, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of substrate/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or egg contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5. The samples should be diluted so that the results can be read in the linear part of the calibration curve.

Please note:

When working according to the described sample preparation, the dilution factor is 20. The whole egg powder concentration can be read directly from the standard curve (the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. *)).

Sample dilutions of more than 1:20 must be considered for the calculation of the whole egg powder concentration.

The test is calibrated against NIST Reference Material 8445 (whole egg powder). The result is expressed as mg/kg whole egg powder. The reference material 8445

contains 49 % total protein. If the result is multiplied by 0.49, then the result can be expressed as mg/kg (ppm) total protein. If the result is multiplied by 0.263 then mg/kg (ppm) egg white protein are obtained.

Calculation example

From the standard curve 10 mg/kg whole egg powder are read. If multiplied by 0.49 then 4.9 mg/kg (ppm) whole egg protein is obtained. If multiplied by 0.263 then 2.63 mg/kg (ppm) egg white protein is obtained.

Recommendation:

In order to ensure a high analytical performance:

- each sample material should be analyzed in duplicates
- egg free and egg containing (spiked) samples should be run as test controls

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG