Peak Flow Meter

BS EN ISO 23747:2007



Vitalograph®

Medical Devices Directive
93/42/EEC L169, Vol. 36.
EN ISO 13485
BSI MD82182 & CE 00772
FDA QSR 21 CFR 820/803

Manufacturer: Vitalograph (Ireland) Ltd, Ennis, Ireland

English

Warnings

Note: Please read all the information in this manual before using the Vitalograph Peak Flow Meter.

- The Vitalograph Peak Flow Meter is recommended for single patient use.

- The Vitalograph Peak Flow Meter should be used under the supervision of a licensed healthcare professional.
- Note: Patient symptoms take precedence over Peak Flow Meter readings.
- If the Vitalograph PFM is used for longer than its specified life, the accuracy of the device may deteriorate.

What is a peak flow meter?

A peak flow meter is a device intended to measure how well your asthma is under control. The device measures air flowing out of the lungs as you blow into the device as fast as possible. A peak flow meter can reveal narrowing of the airways well in advance of an asthma attack. Used mainly by persons with moderate to severe and persistent asthma, peak flow meters can help determine:

- when to seek emergency medical care.
- the effectiveness of a person's asthma management and treatment plan.
- when to stop or add medication, as directed by your physician.
- what triggers the asthma attack (such as exercise-induced asthma).

With asthma, sometimes you may feel your breathing is fine, but when you measure it with a peak flow meter, your lung function may decrease. A peak flow meter can help you determine airways change and better manage your asthma.

Asthma Care in Partnership with your Physician

Your doctor will take the time to explain to you the self-management of your asthma. This will start upon diagnosis and continue with all members of the health team. Your Action Plan will be tailored to your needs, but will include basic facts about asthma, rules of medication, skills required for you to inhaler spacer and Peak Flow Meter, environmental control measures, how to take medicine and how to use a bronchodilator.

With asthma, sometimes you may feel your breathing is fine, but when you measure it with a peak flow meter, your lung function may decrease. A peak flow meter can help you determine airways change and better manage your asthma.

You and Your Asthma

The most people with asthma need to monitor their asthma at least twice a day and to have a plan of action to keep it under control. This leaflet provides detailed information on your Peak Flow Meter and how the Self-Management Action Plan works.

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Réponse: Un Action Plan est une mesure de how fast you can blow air out of your lungs. Your peak flow score will be higher when you are well, and lower when your airways get narrower. Peak flow scores are therefore a useful guide to the openness of your airways at any given moment and hence, the severity of your asthma. You should measure peak flow as soon as you wake up and in the evening before taking your bronchodilator, or as directed by your doctor.

Although how you feel and what you do is important, peak flow scores show accurately how your breathing is changing. Modern asthma medicines aim to give you the best possible peak flow score, keeping your airways open. It is also important to aim for stable peak flow scores over several days. Your treatment and/or Action Plan may be changed following the diagnostic phase. This procedure may be repeated until your optimum Action Plan is proven.

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Q. How do I get an Action Plan?

Réponse: Only your doctor can decide the best Action Plan for you. This is likely to be preceded by an initial assessment followed by a diagnostic phase. During the diagnostic phase you will need to review your peak flow scores [3]. Your Action Plan is then assessed against your peak flow scores over several days. Your treatment and/or Action Plan may be changed following the diagnostic phase. This procedure may be repeated until your optimum Action Plan is proven.

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You should visit your doctor at least twice a year to assess any changes in your condition and change your plan if required. If you find yourself in the yellow zone most of the time, inform your doctor immediately.

Q. Your Action Plan

Self-management including peak flow monitoring will significantly improve the outcome of your asthma.

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Q. Your measured best Peak flow is _____ l/min BTSP

Your Best means the highest Peak flow when you are feeling well.

Q. Updating Your Action Plan

You should visit your doctor at least twice a year to assess any changes in your condition and change your plan if required. If you find yourself in the yellow zone most of the time, inform your doctor immediately.

Q. Understanding your Action Plan

Your current Action Plan is written down and on a label (if provided) which is attached to the back of your Peak Meter. Your doctor or specialist nurse may also mark your "zones" on the scale side of your peak flow meter. If so, use these "colour bands" instead of the scale reading in your day to day management.

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