RIDASCREEN[®] FAST Lysozym

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RIDASCREEN[®]FAST Lysozym

Brief information

RIDASCREEN[®]FAST Lysozym (Art. No.: R6452) is a sandwich enzyme immunoassay for the quantitative analysis of lysozyme (hen's egg protein) in food like wine, cheese and sausage.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation..... approx. 15 min
test implementation (incubation time)..... 30 min

Limit of detection: Wine 0.02 mg/kg (ppm) lysozyme
Cheese, sausage..... 0.10 mg/kg (ppm) lysozyme

Limit of quantification: Wine 0.05 mg/kg (ppm) lysozyme
Cheese, sausage..... 0.25 mg/kg (ppm) lysozyme

Specificity: The antibodies specifically detect lysozyme of hen's egg.

1. Intended use

RIDASCREEN®FAST Lysozym (Art. No.: R6452) is a sandwich enzyme immunoassay for the quantitative analysis of lysozyme (hen's egg protein) in food like wine, cheese and sausage.

2. General

Lysozyme is due to its bacteriostatic effect used for the preservation of food products (E1105). It is one of the allergenic egg white proteins and is contained to 3.5 % in egg white protein. In addition to ovalbumin, ovomucoid and ovotransferrin lysozyme can lead people to allergic reactions. Under the European Allergen Regulation egg has to be declared in food.

Lysozyme is used as a preservative in wine production as well as in cheese and sausage production.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against lysozyme. By adding the standard or sample solution to the wells, present lysozyme will bind to the specific antibodies. The result is an antibody-antigen-complex.

In a washing step components not bound are removed. Then antibody conjugated to peroxidase (enzyme conjugate) is added. This antibody conjugate is bound to the antibody-antigen-complex. An antibody-antigen-antibody-complex (sandwich) is formed.

Substrate/chromogen is added after removing of any unbound enzyme conjugate in a washing step. Bound enzyme conjugate converts the chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is proportional to the lysozyme concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with anti-lysozyme antibodies
- 5 x Standards ^{*)}, 1.3 ml each
0 ppm (zero standard), 0.05 ppm, 0.1 ppm, 0.2 ppm, 0.4 ppm lysozyme
in aqueous solution, ready to use
- 1 x Conjugate (7 ml)red cap
peroxidase conjugated antibody, ready to use
- 1 x Red Chromogen Pro (10 ml) brown cap
substrate/chromogen, stained red
- 1 x Stop solution (14 ml)yellow cap
contains 1 N sulfuric acid
- 1 x Allergen Extraction buffer (125 ml)green cap
as a **10fold concentrate**
- 1 x Wash buffer (100 ml) brown cap
as a 10fold concentrate

*) **The dilution factor 20** which results after sample preparation of wine has already been considered for the standard concentrations. Therefore, the lysozyme concentrations of wine samples can directly be read from the standard curve.

In the case of cheese and sausage the **dilution factor 100** results. Therefore, the concentration read from the standard curve has to be multiplied additionally by factor 5.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifugal vials
- water bath (60 °C / 140 °F)
- fluted paper filter
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The washing buffer contains thimerosal. Avoid skin contact.

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of samples

9.1. Preparations

Improper laboratory equipment can result in contamination of samples or test components with lysozyme. Therefore, gloves should be worn during the test procedure, and before the start of the test note the following:

- thoroughly clean surfaces, glassware, grinders and other equipment
- carry out sample preparation and ELISA assay in separate rooms

The **Allergen Extraction buffer** is available as a **10x concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate **1:10 (1+9) with distilled water** before use (e.g. 100 ml of buffer concentrate + 900 ml distilled water, or alternatively fill up the entire contents of the bottle with dist. water to 1250 ml). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. twelve weeks.

For the extraction of cheese and sausage samples, a high-salt allergen extraction buffer is needed separately. For this 10 g NaCl is added to 100 ml Allergen Extraction buffer. The buffer can be used for 4 weeks.

A representative amount of sample (5 - 50 g) should be homogenized.

9.2. Sample extraction

9.2.1. Cheese and sausage

- add 1 g of a representative, homogeneous sample to 20 ml high-salt allergen extraction buffer (preheated to 60 °C / 140 °F) and mix
- incubate samples for 10 min at 60 °C (140 °F) in a water bath
- cool samples to room temperature (20 - 25 °C / 68 - 77 °F) (e.g. in an ice bath)
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39.2 °F) and/or filter the extract
(alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute supernatant 1:5 (1+4) with finally diluted Allergen Extraction buffer (no salt added)
(**Note:** use this diluted supernatant directly in the test or the diluted extract should be used within 2 hours)
- use 100 µl per well in the test

If further dilutions are necessary, they should be performed with a mixture of Allergen Extraction buffer / high-salt allergen extraction buffer in proportion 5:1 (4 ml + 1 ml).

9.2.2. Wine

- add 1 ml of a representative, homogeneous sample to 19 ml of Allergen Extraction buffer (preheated to 60 °C / 140 °F) and mix
- incubate samples for 10 min at 60 °C (140 °F) in a water bath
- cool samples to room temperature (20 - 25 °C / 68 - 77 °F) (e.g. in an ice bath)
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39.2 °F) and/or filter the extract
(alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl per well in the test

When using wine with high **polyphenol content** (e.g. red wine) add 0.25 % casein (e.g. SIGMA, C5890) or 0.25 % gelatin (e.g. SIGMA, G-6144) to the diluted Allergen Extraction buffer. Add e.g. 0.25 g of casein or gelatin to 100 ml final diluted Allergen Extraction buffer and stir for 10 minutes. Add 19 ml of diluted allergen extraction buffer with additive to 1 ml of wine and extract for 10 min at 60°C.

The extract can be stored 2 days at 4 °C (39.2 °F) or in small portions several months at - 20 °C (- 4 °F).

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Washing Buffer** is provided as a 10fold concentrate. Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.

2. Add 100 µl of each standard solution or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted washing buffer (see 10.1.) and pour out the liquid again. Repeat three more times.
4. Add 100 µl of the enzyme conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted washing buffer (see 10.1.) and pour out the liquid again. Repeat three more times.
6. Add 100 µl of the reddish Substrate/Chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or lysozyme contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5. The samples should be diluted so that the results can be read in the linear part of the calibration curve.

Please note:

When working in accordance with the regulation stated for wine, the dilution factor is 20. The lysozyme concentration can be read directly from the standard curve -

the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. *).

For sample dilutions of more than 1:20 (sample preparation of cheese and sausage) the further dilution factor must be considered for the calculation of the lysozyme concentration.

In case of cheese and sausages the additional factor is 5.

Recommendation:

In order to ensure a high analytical performance:

- each sample material should be analyzed in duplicates
- use also lysozyme free and lysozyme containing (spiked) samples as test controls

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG