RIDA[®] QUICK Zearalenon RQS

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RIDA[®]QUICK Zearalenon RQS

Brief information

RIDA[®]QUICK Zearalenon RQS (Art. No.: R5504) is an immunochromatographic test for detection of zearalenone in corn.

All reagents required for the assay are contained in the test kit.

Results are evaluated with RIDA[®]QUICK SCAN (RQS, Art. No.: ZG5005).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)..... approx. 10 min
test implementation (incubation time)..... 5 min

Limit of detection approx. 75 µg/kg (ppb)

Specificity: The RIDA[®]QUICK Zearalenon RQS test detects zearalenone in corn samples.
The test has been validated for this matrix.

1. Intended use

RIDA[®]QUICK Zearalenon RQS is a quantitative (evaluation with RIDA[®]QUICK SCAN) immunochromatographic test in strip format for the detection of zearalenone in corn.

2. General

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*.

Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hyperestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cow, horse and sheep, too.

The potential health risk for man induced by this mycotoxin, which is taken up with foods of vegetable or animal origin, is extensively discussed.

3. Test principle

The basis of the immunochromatographic assay in test strip format is an antigen-antibody reaction. A specific antibody against zearalenone binds the zearalenone molecules in the samples.

Results are evaluated with RIDA[®]QUICK SCAN (RQS) reader.

The control band (control line) is not influenced by zearalenone in the sample and should be present in all cases in order to prove the test strip is valid.

4. Reagents provided

Each kit contains sufficient materials for 20 determinations.

Each test kit contains:

20 x test strips (one for each determination, separately packed)

1 x mobile solvent (13 ml)

5. Materials required but not provided

5.1. Equipment:

- laboratory mincer / grinder
- balance
- optional: shaker
- optional: centrifuge
- graduated cylinder
- 100 µl, 600 µl micropipettes (e.g. R-Biopharm FP 100, Art. No.: Z0007)
- RIDA[®]QUICK SCAN (RQS, Art. No.: ZG5005)

5.2. Reagents:

- methanol (70 %, (v/v))
- distilled water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively, therefore keep the strips away from humidity. This is especially important for already opened test strip packing.

7. Storage instructions

Store the kit at 2 - 8 °C / 36 - 47 °F. DO NOT FREEZE the test strips.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Allow test strips and buffer to reach room temperature (20 - 25 °C / 68 - 77 °F) before use!

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1. Detection limit 75 µg/kg (ppb)

- weigh 10 g of ground sample into a screw cap tube and add 20 ml of methanol/water (70/30) *)
- close the tube and shake the sample vigorously for approx. 5 min (manually or with shaker)
- let the solution come to sedimentation for approx. 3 - 5 min or centrifuge
- mix 100 µl of the clear supernatant with 600 µl mobile solvent
- apply 100 µl of diluted sample onto the application area of the test strip

*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly, e.g. 50 g in 100 ml of methanol/water (70/30)

9. Test procedure

- apply 100 µl of the sample extract onto the application area of the test strip
- read the result after 5 min (+ max. 15 sec) incubation with RIDA[®]QUICK SCAN

10. Evaluation with RIDA®QUICK SCAN (ZG5005)

Please first read the guide for RIDA®QUICK SCAN attentively.


Please find a short description regarding the use of RIDA®QUICK SCAN in the test kit.

The right band in the field of reaction is a control band (control line) and must appear after each test procedure. If the band is missing, the test result is not valid because of improper test procedure or deterioration of the reagents. Repeat the test with a new strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is contaminated > 1000 ppb the control band is visible very weak.

To obtain comparable results, the test strips must be evaluated always after 5 min.

Fig.1: test strip RIDA®QUICK Zearalenon RQS

11. Sensitivity

The RIDA®QUICK Zearalenon RQS test is capable of detecting zearalenon contaminations of $\geq 75 \mu\text{g/kg}$ (ppb) to $500 \mu\text{g/kg}$ (ppb).

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG