# **RIDASCREEN<sup>®</sup> FAST Soya**

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# RIDASCREEN<sup>®</sup>FAST Soya

## Brief information

RIDASCREEN<sup>®</sup>FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of soya proteins in untreated and processed food and beverages.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) ..... approx. 25 min  
test implementation (incubation time) ..... 30 min

Kalibrator: soya protein (soya flour with 39 % soya protein)

Limit of detection: 0.31 mg/kg (ppm) soya protein

Limit of quantification: 2.5 mg/kg (ppm) soya protein

Specificity: The antibodies specifically detect heated soya proteins.  
(The sample extraction at **100 °C for 10 min** ensures that all soya proteins denature and can be detected in the ELISA).

There is a weak cross-reactivity of 0.0017 % to legumes of the tribe Phaseoleae (various species of beans) and to legumes of the genus *Vicia* (0.0003 %).

There is no cross-reactivity to other legumes like peanut, lentil, pea, lupine as well as milk (casein,  $\beta$ -lactoglobulin) or egg proteins.

## 1. Intended use

RIDASCREEN®FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed soya proteins in food. Such products are beverages as juice, wine, beer and food as sausages, dressings, bakery products, ice, chocolate, soup, sauce, margarine, etc.

## 2. General

Soya proteins can be present as an ingredient or as a contamination in raw and processed products. According to the EU Directive 2003/89/EG from November 10<sup>th</sup> 2003, soya and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Soybeans (mature seeds, raw) contain approx. 35 % proteins. Therefore, soya is one of the few plants which are considered as a valuable source of proteins. Thus, soya is widely used as substitute for proteins from animal sources (e.g. in tofu, soya milk or soya yoghurt). Besides the high protein content, soya beans are also rich in fat (approx. 20 %) and are used for the productions of oils and fats.

The usage of soya as food and feed has strongly increased in the last decades, leading to an increase in allergic reactions as well.

## 3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against soya proteins. By adding the standard or sample solution to the wells, present soya proteins will bind to the specific antibodies. The result is an antibody-antigen-complex.

In a washing step components not bound are removed. Then antibody conjugated to peroxidase (enzyme conjugate) is added. This antibody conjugate is bound to the antibody-antigen-complex. An antibody-antigen-antibody-complex (sandwich) is formed. Substrate/Chromogen solution is added after removing of any unbound enzyme conjugate in a washing step. Bound enzyme conjugate converts the chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is proportional to the soya protein concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)  
coated with anti-soya protein antibodies
- 5 x Soya protein standard solutions <sup>\*)</sup>, 1.3 ml each  
0 mg/kg (ppm) (zero standard), 2.5 mg/kg, 5 mg/kg, 10 mg/kg, 20 mg/kg  
soya protein in aqueous solution, ready to use
- 1 x Conjugate (0.7 ml) .....red cap  
peroxidase conjugated antibody, 11x concentrate
- 1 x Red Chromogen Pro (10 ml) ..... brown cap  
Substrate/Chromogen solution, stained red
- 1 x Stop solution (14 ml) ..... yellow cap  
contains 1 N sulfuric acid
- 1 x Allergen Extraction buffer (125 ml) .....green cap  
as a 10x concentrate
- 1 x Washing buffer (100 ml) ..... brown cap  
as a 10x concentrate
- 4 x Extractor 3 (30 ml) ..... brown cap  
Extraction solution for samples

\*) The dilution factor of 100 for the sample has already been considered when labeling. Therefore, the soya protein concentration of samples can directly be read from the standard curve.

## 5. Reagents required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge, centrifugal vials
- shaker
- water bath (60 °C / 140 °F and 100 °C / 212 °F)
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- paper filter
- graduated pipette
- variable 20 - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistepper pipette

## 5.2. Reagents:

- distilled or deionized water
- for the sample preparation of tannin and polyphenol containing food: casein (e.g. Sigma C5890) and polyvinylpyrrolidon (PVP) K15 or K17 (e.g. AppliChem A2258)

## 6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instructions for use must be strictly followed.

Extractor 3 is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

The washing buffer contains thimerosal. Avoid skin contact.

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

Red Chromogen Pro is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of Red Chromogen Pro prior to test implementation
- a value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for standard 5
- odorless Extractor 3 after repeated opening of the bottle; Extractor 3 should smell like mercaptoethanol (rotten eggs). (For checking, carefully fan the vapor of the open bottle)

## 9. Preparation of Samples

The **Allergen Extraction Buffer** is provided as a **10fold concentrate**. Before dilution, dissolve possibly formed crystals by heating (water bath 37 °C / 98.6 °F) and mix well. Dilute the heated concentrate 1:10 (1+9) with distilled water (e.g. 100 ml concentrate + 900 ml dist. water) or fill up the whole content of the bottle to 1250 ml. The diluted Allergen Extraction Buffer can be used for approx. 12 weeks at 2 - 8 °C (35.6 – 46.4 °F).

### 9.1. Sample extraction

Unclean laboratory equipment can result in a contamination of samples or test components with soya. Therefore, gloves should be worn during the test procedure and the following recommendations should be applied:

- thoroughly clean surfaces, glass vials, mixers and additional equipment before and after each sample
- sample extraction and ELISA-procedure should be performed in separated rooms
- homogenize a representative amount of the sample (5 - 50 g)

### 9.2. Solid samples:

- weigh 1 g of the homogenized sample, add 2.5 ml Extractor 3 and 17.5 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction Buffer, close the vial, mix vigorously and put the solution for **10 min at 100 °C / 212 °F** in a boiling water bath
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with diluted Allergen Extraction Buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction Buffer)
- use 100 µl extract per well in the assay

### 9.3. Liquid samples

- take 1 ml of the homogenized sample, add 2.5 ml Extractor 3 and 16.5 ml preheated (60 °C / 140 °F) and diluted Allergen Extraction Buffer, close the vial, mix vigorously and put the solution for **10 min at 100 °C / 212 °F** in a boiling water bath
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract

(alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)

- dilute the extract 1:5 (1+4) with diluted Allergen Extraction Buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction Buffer)
- use 100 µl extract per well in the assay

#### 9.4. Tannin and polyphenol containing food (spices, chocolate, coffee)

To the finally diluted Allergen Extraction Buffer (AEB), add 0.25 % casein and 1 % PVP, this aids extracting soy proteins from different spices, chocolate and other phenol containing matrices (see validation data, available at R-Biopharm). Components are better soluble if the solution is preliminarily heated to 60°C / 140 °F.

- weigh 1 g of the homogenized sample, add 2.5 ml Extractor 3 and 17.5 ml preheated (60 °C / 140 °F), finally diluted AEB (containing 0.25 % casein and 1 % PVP), close the vial, mix vigorously and put it for **10 min at 100 °C / 212 °F** in a boiling water bath
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with finally diluted Allergen Extraction Buffer (e.g. 100 µl sample extract + 400 µl finally diluted Allergen Extraction Buffer)  
use 100 µl extract per well in the assay

**The diluted Casein/PVP-Extraction Buffer can be stored at 2 - 8 °C (35 - 46 °F) for 7 days. The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. The extracts can be stored at -20 °C (-4 °F) for several months.**

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **antibody enzyme conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted enzyme conjugate solution has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **washing buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for four weeks.

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at the same time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard solution or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted enzyme conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 -77 °F).

5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl Red Chromogen Pro to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ( $A_{450\text{ nm}}$ ) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or soya contamination.

A further dilution and new detection of the samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 5. The samples should be diluted so that the results can be read in the linear part of the calibration curve.

### **Please note:**

When working in accordance with the regulation stated, the dilution factor is 100. The soya protein concentration can be read directly from the standard curve - the sample dilution factor of 100 is already taken into account for the standard concentrations (see 4. \*).

For sample dilutions of more than 1:100 the further dilution factor must be considered for the calculation of the soya protein concentration. E.g. if a sample is tested for 1 % soya flour concentration (corresponding to 10000 mg/kg soya flour or 3900 mg/kg soya protein), the sample has to be diluted 1:70000 (1:100 according to the sample preparation and additional 1:700 with Allergen Extraction Buffer).

**Remarks:**

A result of 10 mg/kg soya protein corresponds to 25.6 mg/kg soya flour (in case of 39 % protein content).

Samples tested negative still could contain:

- a soya contamination below the limit of detection of the assay
- protein fragments (hydrolysates, isolates) which are not fully determined in the sandwich format.

**Recommendation:**

In order to ensure a high analytical performance:

- each sample material should be measured in duplicates
- employ also soya protein free and soya protein containing (spiked) samples as test controls

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG