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Set of 3 processed Gliadin Assay Controls

1. Intended Use

The processed assay controls are used to control the extraction, the test implementation and the handling of the RIDASCREEN® Gliadin (R7001) and RIDASCREEN®FAST Gliadin (R7002).

It is recommended to prepare the assay controls with Cocktail (patented) (Art. No. R7006/R7016, official R5-Mendez method).

2. Content

Each assay control contains 1.5 g homogenized snack samples.

Gliadin Assay Control A:	approx. 5.5 ppm prolamin
Gliadin Assay Control B:	approx. 20 ppm prolamin
Gliadin Assay Control C:	approx. 50 ppm prolamin

The Set of 3 processed assay controls is composed of 3 snack samples with a defined content of gliadin. These assay controls were produced with a highly processed and ground snack, exactly characterized in consideration of its gliadin-content^[1]. Furthermore, these defined snack samples were analyzed in an official interlaboratory study of AACCI and PWG. For the preparation of assay controls, these snack samples were partly remixed.

3. Storage

Store sample at 2 - 8 °C (35.6 - 46.4 °F) in the dark.

4. Preparation of samples

It is recommended to prepare the samples with Cocktail (patented, R7006/R7016), as explained in the product information of RIDASCREEN® Gliadin (R7001) and RIDASCREEN®FAST Gliadin (R7002).

5. Test implementation of the ELISA

The further procedure is explained in the product information of RIDASCREEN® Gliadin (R7001) and RIDASCREEN®FAST Gliadin (R7002).

6. Results

The characteristic of RIDASCREEN® Gliadin and RIDASCREEN®FAST Gliadin is the quantitative analysis of prolamins in raw products as well as in processed food.

A gliadin concentration in µg/kg (ppb) is read out from the RIDA®SOFT Win calibration curve and must be multiplied with the dilution factor of 500 or 2000.

If the results do not correlate with the values in the table, a trouble shooting of the ELISA procedure has to be carried out.

Often airborne cereal dust or dirty laboratory equipment lead to gliadin contamination in the assay. It is recommended to clean surfaces, pipettes, glass vials and mincers before use with 40 % ethanol or propanol.

Assay control / Test Kit *	Dilution factor (final)	MV * (ppm Gliadin)	CV * (%)	Range (ppm Gliadin)
A				3.8 – 10.0
R7001 (n=10; r=30)	500	6.3	15.8	
R7002 (n=10; r=30)	500	5.5	16.5	
B				16.5 – 30.6
R7001 (n=10; r=30)	500	23.8	10.2	
R7002 (n=10; r=30)	500	22.1	10.9	
C				37.0 – 68.6
R7001 (n=10; r=30)	2000	53.2	8.2	
R7002 (n=10; r=30)	2000	48.9	7.4	

Assay Control / Test Kit *: n = number of extractions; r = number of repeated measurements

MV *: mean value

CV *: coefficient of variation

Remark

The gliadin assay controls are **not** a reference material. For validation of the RIDASCREEN® Gliadin tests spiked samples should be used, applications can be provided by R-Biopharm.

Literature

[1] Koehler P, Schwalb T, Immer U, *et al.* AACCI Approved Methods Technical Committee Report: Collaborative study on the immunochemical determination of intact gluten using an R5 sandwich ELISA. *Cereal Foods World* 2013; **58**: 36-40.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG