



# Papanicolaou OG6

## Cytoplasmic staining



In-vitro diagnostic medical device

**BASIC UDI: 080339762W01030708X8**

IVD in **Classe A**, Reg. UE 2017/746

Catalog number	Unit size	UDI-DI
05-12013	500 ml	08033976232702
05-12013/L	1 l	08033976232719
05-12013E	2,5 l	08033976232726

### Packaging

- 05-12013E

Primary container: white bottle in polyethylene terephthalate (PET). Useful capacity 2.5 liters. HDPE cap.

Tamper evident cap.

The polyethyleneterephthalate is a thermoplastic polymer of the polyester family. PET is an optimal oxygen, carbon dioxide and other gasses barrier. This material has a high resistance to ultraviolet radiation and an inertia toward the mainly chemical agents (solvents: xylene, limonene, liquid paraffines, alcohols, acids, bases etc.). It is biologically inert. It constitutes a good water and humidity barrier. It shows a great hardness and mechanical resistance.

The bottle has an optimal grip. The absence of the handles reduces space for storage. The anti-dropping cap permits a precise and clean use.

Secondary container: carton box.

- 05-12013/L

Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 1 l. HDPE cap. Tamper evident cap.

-05-12013

Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 500 ml. HDPE cap. Tamper evident cap.

### Expected aim

Product for the preparation of: gynecological specimens, urine cytology, fine needle specimens, sputum and bronchial washings, to be examined by optical microscopy.

### Application

Cytoplasmic staining for the keratinized cells in the Papanicolaou method.

**For the execution of the staining method is required the use of reagents Papanicolaou Harris hematoxylin and Papanicolaou EA50.**

**Principle** A highly selective blue nuclear stain, Harris' hematoxylin, is combined with EA50 polychromic mixture, a subtle cytoplasmic stain which differentiates cyanophil cells from eosinophil ones. The last ingredient is OG6 solution, which stains keratinized elements.

**Method**

- 1) Ethanol 95°, 2 minutes
- 2) Distilled water, 2 minutes
- 3) Harris Hematoxylin, 1 minute
- 4) Tap water, 5 minutes
- 5) Ethanol 95°, 15 seconds
- 6) OG 6, 2 minutes
- 7) Ethanol 95°, 15 seconds (twice)
- 8) EA 65, 5 minutes
- 9) Ethanol 95°, 15 seconds
- 10) Absolute Ethanol, 30 seconds (twice)
- 11) Xilene or Bio Clear, 2 minutes (twice)

**Results**

Nuclei : Blue-purple  
 Cyanophil cytoplasm : Blue-green  
 Eosinophil cytoplasm : Pink  
 Keratinized cytoplasm : From pink to orange

Components	Components	CAS	CE	Index
	Orange G	1936-15-18	-	-
	Phosphotungstic acid	12501-23-4	-	-
	Ethanol 95°	64-17-5	200-578-5	603-002-00-5
	Deionized water			

**Warning and precaution**

The product is intended for professional laboratory use for healthcare professionals. Carefully read the information on the label (danger symbols, risk and safety phrases) and always consult the safety data sheet. Do not use if the primary container is damaged. In the event of a serious accident, we recommended that you immediately inform Bio-Optica Milano S.p.A and the competent authorities.

**Storage** Store the preparation at 15-30°C. Keep the containers tightly closed.

**Stability** After the first opening, the product is reusable until the expiry date, if correctly stored. Validity: 2 years.

**Disposal** Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

**References**

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022