

## Immersion oil for microscopy

**IVD** In-vitro diagnostic medical device **CE**

**CND Code: W01030704**

Catalog number	Unit size
08-1730/A30	30 ml
08-1730/A270	9 x 30 ml

**Packaging** Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 30 ml. Precision dropper in HDPE. Tamper evident cap in HDPE.  
Secondary container: white carton box + label.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

**Expected aim** Product for the preparation of cyto-histological samples for optical microscopy.

**Specifications** Refractive index: 1,5  
Density: 862 Kg/m<sup>3</sup> (20°C)  
Viscosity: 60 cSt (40°C)

**Application** Oil type A for optical microscopy

**Method** Put one drop on the smear and then focus by the immersion objective.  
To clean the objectives it is recommended to use of the Bio-Optica cloth (cod. 09-000500)

**Components**

Components	CAS	CE	Index
Mixture of paraffinic oils	-	-	-

**Functional properties** The product dries slowly and this enables to clean the objective also after many days after its use.

**Warning and precaution** The product must be used exclusively by specialized technical operators.  
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.  
Do not use if the primary container is damaged.

**Storage** Store in a cool dry area, at room temperature. Keep away from sources of ignition.  
Keep from light . Keep the containers tightly closed.

**Stability** After the first opening, the product is usable until the expiry date, if correctly stored. Validity: 2 years.

**Disposal** Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: July 2013