# **RIDA<sup>®</sup> QUICK Gliadin**

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# RIDA<sup>®</sup> QUICK Gliadin

## Brief information

RIDA<sup>®</sup>QUICK Gliadin (Art. No.: R7003) is an immunochromatographic test for the qualitative detection of Gliadin/Gluten contamination on surfaces (swab test for the hygiene control in the production and in laboratories) as well as in gluten-free raw material.

All reagents required for the assay are contained in the test kit.

The test kit is sufficient for 25 determinations.

Results are read visually.

Time requirement:      sampling for swab test..... approx. 1 min  
                                 sample preparation for 10 raw materials ... approx. 15 min  
                                 (homogenization, extraction, centrifugation)

                                 test implementation (incubation time)..... 5 min

Detection limit:      on **surfaces** approx. 0.5 µg gliadin / 100 cm<sup>2</sup>  
  
                                 in **raw material** approx. 2.5 mg/kg gliadin (approx. 5 mg/kg  
                                 gluten)

Specificity:              The monoclonal antibody R5 reacts with the gliadin-fraction  
                                 from wheat and corresponding prolamins from rye and  
                                 barley.  
                                 No cross-reaction with soy, oats, corn, rice, millet, teff,  
                                 buckwheat, quinoa and amaranth.

## Related products:

RIDASCREEN<sup>®</sup> Gliadin (Art. Nr.: R7001)

RIDASCREEN<sup>®</sup>FAST Gliadin (Art. Nr.: R7002)

RIDASCREEN<sup>®</sup> Gliadin competitive (Art. Nr.: R7011)

RIDA<sup>®</sup> Extraction Solution (Art. Nr.: R7099)

Cocktail Solution (Art. Nr.: R7006)

Set of 3 Gliadin Assay Controls (Art. Nr.: R7010)

SureFood<sup>®</sup> Allergen real time PCR Gluten (Art. Nr.: S3106)

## **1. Intended use**

RIDA<sup>®</sup> QUICK Gliadin can be used as a swab test for the gluten determination on surfaces in the hygiene control and for the qualitative detection of gliadin/gluten in raw material. The test should only be used for the detection of small amounts of gluten (contaminations).

The gliadin/gluten determination can be carried out for numerous raw materials after a simple ethanol extraction. Further applications, e.g. for processed foods and chocolate, are available upon request.

## **2. General**

The use of wheat flour and gluten in foodstuff is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

In the document of Codex Alimentarius (Alinorm 08/31/26) gluten-free foods are defined as foods, which having a gluten content of less than 20 mg/kg.

## **3. Test principle**

The immunochromatographic test employs the monoclonal R5-antibody which is specific for the detection of gliadin from wheat and prolamins from rye and barley. Results are read visually. Generally, the higher the analyte level in the sample, the stronger the red color of the test band will be.

## **4. Reagents provided**

Each kit contains sufficient materials for 25 determinations. Each test kit contains:

25 x dip sticks (one for each determination)

25 x disposable pipettes

30 x test tubes

1 x sample dilution buffer (60 ml), 5x concentrate

## 5. Materials required but not provided

### 5.1. Equipment:

#### **For swab tests**

- graduated cylinder for preparing the sample buffer
- vessel for preparing the sample dilution buffer

#### **For analysis of raw material**

- balance
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- shaker
- centrifuge + centrifugal glass vials
- graduated pipettes

### 5.2. Reagents:

#### **For swab tests**

- distilled or deionized water

#### **For analysis of raw material**

- distilled or deionized water
- ethanol solution (60 %), for the extraction of the samples (add 150 ml ethanol p.a. to 100 ml distilled water and shake well)
- for soy containing food: skim milk powder (food quality)

## 6. Warnings and precautions for users

Airborne cereal dust and dirty laboratory equipment lead to gliadin contamination of the assay. In order to avoid cross-contamination by cereal dust, please note the following points:

- wear gloves before starting and during the assay
- clean surfaces, glass vials, mincers and other equipment with 60 % ethanol
- for raw material the extraction and the test procedure should be carried out in separate rooms

The dip sticks are very sensitive to humidity, that could turn the test useless. For this reason keep the strips away from humidity!

## 7. Storage instructions

Store the unopened kit at 2 - 8 °C (36 - 46 °F). DO NOT FREEZE.

Once the seal of the test strip container is broken, store the container at room temperature (20 - 25 °C / 68 - 77 °F).

No quality guarantee is accepted after the expiration date on the kit label.

## 8. Test implementation

### 8.1. Test preparation

The **sample dilution buffer** is provided as a 5fold concentrate. Only the amount which actually is needed should be diluted 1:5 (1+4) with distilled water (e.g. 1.2 ml concentrate + 4.8 ml distilled water, sufficient for 10 samples). Make sure that the buffer is not contaminated with gliadin from the sample.

### 8.2. Swab test: sampling and test implementation

1. Take as many test tubes as samples to be analysed.
2. Place 500 µl of diluted sample dilution buffer in the test tube (e.g. using the disposable pipette provided).
3. Swab the lower end (reaction zone) of a dry dip stick thoroughly over a sampling area of 10 x 10 cm (wear gloves).
4. Place the dip stick inside vertically with the arrow end into the tube. Do not immerse the strip past the maximum line.
5. Take out the strip after exactly 5 min and read the result immediately without further manipulation.

### 8.3. Sample preparation raw material

#### 8.3.1. Fluid and soft food

- fluid food:** mix 1 ml of the sample with 9 ml 60 % ethanol solution
- for soy milk add additionally 1 g of skim milk powder
- soft food:** weigh 1 g of a representative sample and add 10 ml 60 % ethanol solution
- for soy milk products add additionally 1 g of skim milk powder
- shake well for at least 30 sec. (vortex)
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F)

#### 8.3.2. Solid and hard food

- weigh 5 g sample and grind it to powder
- use 1 g of this powder and add 10 ml 60 % ethanol solution (for soy containing samples add 1 g of skim milk powder)
- shake well for at least 30 sec. (vortex)
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F)

#### 8.3.3. Samples containing inhomogeneous gliadin content (e.g. meat and sausages)

In these matrices gliadin may not be distributed evenly. Therefore, use more of the sample and the corresponding amount of 60 % ethanol.


- e.g. 20 g of a homogeneous sample and add 200 ml of 60 % ethanol solution
- shake well for at least 30 sec. (vortex)
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F)

#### **Remark:**

**All supernatants obtained after centrifugation can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °C) up to four weeks.**

## 8.4. Test implementation for raw material

1. Take as many test tubes as samples to be analysed.
2. Place 500 µl of diluted sample dilution buffer in the test tube (e.g. using the disposable pipette provided).
3. Pipette 50 µl of the sample supernatant or place 3 drops with the provided disposable pipette, vertically dropped in the test tube and shake slightly.
4. Place the dip stick inside vertically with the arrow end into the tube. Do not immerse the strip past the maximum line.
5. Take out the strip after exactly 5 min and read the result immediately without further manipulation.

C = control band (blue)

T = test band (red)

## 9. Results and Sensitivity

### Positive result: two colored bands

The sample is positive if two colored bands (the blue control band and the red test band) are visible within the result window.

**Swab test:** concentration > 0.5 µg gliadin / 100 cm<sup>2</sup>

**Raw material:** concentration > 2.5 mg/kg gliadin (corresponds to 5 mg/kg gluten)

### **Negative result: only the blue control band**

The sample is negative if no red test band is visible within the result window.

**Swab test:** concentration < 0.5 µg gliadin / 100 cm<sup>2</sup>

**Raw material:** concentration < 2.5 mg/kg gliadin (corresponds to 5 mg/kg gluten)

### **Remark:**

- The test strip has been developed for the detection of gluten contamination.
- When testing samples with a high gluten content the detection can be suppressed (“hook-effect”). If there is a suspicion of a too high concentration, the sample should be diluted at least 1:100.
- The sample extraction with ethanol should only be used for raw material that were surely not heated and not processed.
- The Cocktail Solution or the RIDA<sup>®</sup> Extraction Solution has to be used for processed foods in order to detect also the heat-altered prolamins.
- It is recommended to compare the extraction efficiency of ethanol with the Cocktail Solution (R7006) and the RIDA<sup>®</sup> Extraction Solution (R7099).

### **Invalid result: no colored band**

If no blue or red band is visible within the result window after performing the test, the test is considered invalid. (“Hook effect”: High concentrated samples overload the strip and no band will be visible)

### **Note:**

For documentation the upper part of the dip stick marked with “Gluten” together with the test bands must be cut off.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG