

Small incision, big results



## **VERTEBRIS stenosis**

Full-endoscopic, interlaminar decompression  
for lumbar spinal canal stenosis

# VERTEBRIS *stenosis*

Full-endoscopic Spine Instrumentation

## Contents

<b>Full-endoscopic, interlaminar decompression</b>	
■ Introduction	04
■ Positioning	06
■ Determination of interlaminar access	06
■ Performance of interlaminar access	07
■ Ipsilateral, decompression on one side	08
■ Contralateral decompression in over-the-top technique	10
<b>VERTEBRIS stenosis Instrumentation</b>	
■ Endoscope and accessories VERTEBRIS stenosis	12
■ Access and working instruments VERTEBRIS stenosis	13
■ PowerDrive ART1 – Universal motor system	14
■ COMBIDRIVE EN – High-speed motor system	15
■ Surgimax – Radiofrequency generator	16
■ Consumables and replacement parts	17
Literature	18

# VERTEBRIS stenosis



## Introduction

Degenerative stenosis of the lumbar spinal canal with compression of neural elements arise as a result of bony, disk, capsular or ligament structures. Depending on localization and spread, they can lead to classic symptoms in the lower extremities. Pain in the back tends to be attributed to secondary degenerative phenomena, e.g. segmental instability or deformity. There is no clear correlation between the extent of the stenosis shown by imaging and the clinical symptoms. Apart from spinal disk herniations, lateral and central spinal canal stenoses form the most frequent causes.

A surgical procedure may be necessary after conservative measures have been exhausted or neurological deficits occur. When this is the case, the pathology and symptoms must be taken into account and decompression operations, fusions, or a combination of both procedures must be considered. Today, EBM criteria appear to provide certainty that decompression procedures can improve radicular symptoms and neurogenic claudication. The extent of decompression required from a technical perspective and the conditions under which an additional fusion is necessary have not been definitively described.

Conventional decompression operations on the lumbar spine demonstrate good results. However, consequences and problems associated with these operations are known. Attempts were therefore made right from the start of spine surgery to modify existing operating procedures. Up to the present day, the primary focus continues to be on reducing the invasiveness of surgery and improving the intraoperative view.

Minimally invasive techniques can reduce the trauma and consequences due to the operation. At the same time, visualization and illumination during the operation can be optimized. Appropriate instrument sets for decompression of lumbar spinal canal stenosis were developed on the basis of experiences derived from full-endoscopic operations on spinal disk herniations of the cervical and lumbar spine, offering the possibility of endoscopic bone resection. Since a more extensive bone or ligament resection is frequently necessary here, a large endoscope with a correspondingly large intraendoscopic working channel and larger instruments were necessary. Full-endoscopic, interlaminar access is used routinely, while the transforaminal/extraforaminal access is reserved for specific individual cases.



Lateral and central spinal canal stenosis of the lumbar spine.




A range of endoscopes is available to match different pathologies





Spinal canal decompression with interlaminar access



Intraoperative site after decompression

Today, the instrument sets available permit a full-endoscopic approach under visualization, depending on the indication criteria, which is equivalent to conventional operations. While lateral, stenosis with symptoms on one side can be frequently operated using the basic instrument set, the larger Stenosis System can be used to operate on advanced cases or central stenosis. It is always important to consider whether a stabilizing measure is necessary in addition to decompression.

Sebastian Ruetten, M.D.

Martin Komp, M.D.

**Center for Spine Surgery and Pain Therapy**  
 Head: Priv.-Doz. Dr. med. habil. Sebastian Ruetten



**ST. ELISABETH GRUPPE**   
 KATHOLISCHE KLINIKEN RHEIN-RUHR

**Center for Orthopedics and Traumatology  
 of the St. Elisabeth Group – Catholic Hospitals Rhein-Ruhr**  
**St. Anna Hospital Herne/Marienhospital Herne University Hospital/Marien Hospital Witten**  
 Director: Prof. Dr. med. Georgios Godolias

# VERTEBRIS stenosis


## Full-endoscopic, interlaminar decompression

### Positioning

The patient is placed on an operating table with an X-ray transparent top in the prone position with pelvic and thorax support pillows. A C-arm image intensifier is required during the procedure.




Prone position with pelvic and thorax pillows




### Determination of interlaminar access


Using image intensifier control, the access is determined on the basis of anatomical landmarks in the posterior-anterior beam path and taking account of the pathology. The port must be maximally medial in the interlaminar window in order to permit easier lateral access below the obliquely positioned zygoapophyseal joints.




Marking the entry point on the skin



Entry point should be in a maximally medial position



Access below the zygoapophyseal joints should be possible




Skin incision




## Performance of interlaminar access


After determining the entry point on the skin and performance of the skin incision, the dilator is inserted up to the ligamentum flavum or to the zygoapophyseal joints under posterior-anterior image intensifier control. The subsequent procedure is then performed in the lateral beam path. The working sleeve with oblique opening is pushed over the dilator toward the ligament and the dilator is removed. The endoscope is introduced and the ongoing intervention carried out under continuous visualization and irrigation.



Insertion of the dilator, ...



... operating sleeve ...




... and endoscope

# VERTEBRIS stenosis

## Full-endoscopic, interlaminar decompression

### Ipsilateral, decompression on one side

After the access has been created, the bony structures are exposed. It may be helpful to start decompression at the caudal end of the descending facet. Depending on the pathology, decompression is then commenced with resection of parts of the medial descending facet, the cranial and caudal lamina, and the ligamentum flavum. The extent of decompression generally continues cranially at least until the tip of the ascending facet and caudally to half of the pedicle. The medial portions of the ascending facet and the ligamentum flavum are then resected until sufficient decompression of the neural structures can be clearly seen




It may be helpful to start decompression at the caudal end of the descending facet




The extent of bone resection must generally reach from the tip of the ascending facet to the middle of the caudal pedicle



A range of burrs and bone punches is available for bone resection. They can be introduced through the intraendoscopic working channel



Resection of medial portions of the ascending facet




Removal of protruding annulus and osteophytes



cranially, caudally and laterally. In the case of a central stenosis, the ligamentum flavum generally needs to be resected medially to the midline. Finally, it may be necessary to remove protruding annulus parts and osteophytes in the ventral epidural space. If the patient experiences bilateral symptoms of a lateral stenosis, "over the top" access using the undercutting technique to the opposite side is not carried out. An independent contralateral access is used to retain the median portions of the ligamentum flavum and leave the spinal canal untouched here.



Site after ipsilateral decompression






# VERTEBRIS stenosis


Full-endoscopic, interlaminar decompression

## Contralateral decompression in over-the-top technique

If bilateral symptoms occur with a central stenosis, a unilateral approach is carried out with "over-the-top" access using the undercutting technique to the opposite side. For this purpose, bone in the ventral area of the spinous process is resected until the contralateral side can be accessed dorsally up to the dura of the spinal cord. If possible, the ligamentum flavum is initially left in place to protect the dura and bony decompression is again carried out by laminotomy and partial facetectomy. The ligamentum flavum is then completely resected. Finally, the contralateral recess needs to be extended. The decompression is completed when the dura and the spinal nerves have been clearly decompressed.











In general the sealing caps for endoscope and working sleeve should only be used briefly if bleeding obscures visibility since when operations last a long time and the drainage of fluid is prevented without being noticed, the consequences of volume overload and elevated pressure within the spinal canal and the associated and neighbouring structures should not be ignored. An extended and uninterrupted excessive retraction of the neural structures with the working sleeve in a medial direction must be avoided particularly in cranial areas, or only carried out intermittently, in order to avoid the risk of neurological damage. Experience indicates that as with all new techniques there is generally an enhanced risk of problems occurring during the learning curve.





















Site after over-the-top decompression

# VERTEBRIS stenosis

## Endoscope and accessories

Endoscope	
	<p><b>Discoscope</b>                      PANOVIEW Plus Discoscope 20°, working channel ID 5.6 mm, OD 9.3 x 7.4 mm,                      WL 177 mm (892109205)..... 892109205</p>
	<p><b>Endoscope adapter</b>                      for distance control ..... 892009000</p>
	<p><b>Fiber light cable</b>                      D 3.5 mm, WL 3 m .....806635301</p>
	<p><b>Reprocessing Tray</b>                      for transport and sterilization of Vertebriis Stenosis instruments width 10.4" x length                      21.2" x height 4.6" ..... 88090009</p>
Retainer arm systems	
	<p><b>Retainer arm adapter</b>                      for fixing the endoscope in combination with universal retainer arm (898004717)                      or LEYHLA articulated arm (8766.951) ..... 8840.9722</p>
	<p><b>Universal retainer arm</b>                      max. retaining force 90N, joint radius 420 mm ..... 898004717</p>
	<p><b>Clamp socket</b>                      electrically insulated to 4.5 KV AC, for attaching to the operating table                      (standard rail) ..... 8840.9722</p>
	<p><b>LEYHLA articulated arm</b>                      for attaching to the operating table, fixing in place using 2 articulated arms                      is recommended ..... 8766.951</p>


## Access and working instruments










Access instruments	
	<b>Dilator</b> OD 9.4 mm, cannulated, for working sleeve OD 10.5 mm ..... 892209510
<b>Step-dilator Set comprising:</b> ..... 8922095000	
	<b>Dilator</b> OD 3.9 mm, cannulated, for working sleeve ..... 892209505
	<b>Dilator</b> OD 5.9 mm, for working sleeve or dilator OD 7.0 mm ..... 892209507
	<b>Dilator</b> OD 6.9 mm, for working sleeve or dilator OD 9.5 mm ..... 892209508
	<b>Dilator</b> OD 9.4 mm, for working sleeve OD 10.5 mm ..... 892209515
	<b>Working sleeve</b> OD 10.5 mm , WL 120 mm ..... 892209010
	<b>Irrigation adapter</b> OD 10.5 mm ..... 892209310
Working Instruments	
	<b>Tube sheath punches</b>
	<b>Punch</b> , 2 mm cutting width, OD 5.5 mm, WL 340 mm, TL 490 mm ..... 892409020
	<b>Punch</b> , 3.5 mm cutting width, OD 5.5 mm, WL 340 mm, TL 490 mm ..... 892409035
	<b>Kerrison punches</b>
	<b>Kerrison punch 60°</b> 4.5 mm cutting width, D 5.5 mm, WL 350 mm, TL 460 mm ..... 892409445
	<b>Kerrison punch 90°</b> 4.5 mm cutting width, D 5.5 mm, WL 350 mm, TL 460 mm ..... 892409945
	<b>Micro punches and rongeurs</b> Color coding for easy identification of instrument diameter
	<b>Rongeur</b> OD 3 mm, WL 290 mm ..... 89240.3003
	<b>Rongeur</b> OD 4 mm, WL 290 mm ..... 89240.3004
	<b>Punch</b> OD 3 mm, WL 290 mm ..... 89240.3023
	<b>Punch</b> OD 4 mm, WL 290 mm ..... 89240.3024



# VERTEBRIS stenosis


## PowerDrive ART1 – Universal motor system









Universal motor system	
<b>Burrs for Power Stick M 5</b>	
	<b>Oval burr, with side guard, OD 5.5 mm, WL 290 mm</b> .....899751505
	<b>Oval burr, eccentric, with side guard, OD 5.5 mm, WL 290 mm</b> .....899751555
	<b>Round burr, OD 5.5 mm, WL 290 mm</b> .....899751305
	<b>Round burr, diamond, OD 5.5 mm, WL 290 mm</b> .....899751405
<b>Articulated burr – TipControl</b>	
	<b>TipControl® – Articulating bone burr, complete, OD 4.0 mm, WL 290 mm outer sleeve (899753754) and inner sleeve (15336058) for Power Stick M5, set including burr insert, sterile, pack of 5 (499751704), key for inserting and removing the burr insert (15372005), irrigation adapter (15261106)</b> .....899753794
<b>Motor handles – Power Stick M5</b>	
	<b>Power Stick M5/0</b> Handle for shaver blades or burrs, operation with footswitch, sterilizable, max. speed 16,000 rpm, with fixed connection cable .....89955.0000
	<b>Power Stick M5/3</b> Handle for shaver blades or burrs, operation with keypad or footswitch, sterilizable, max. speed 16000 rpm .....89955.0003
	<b>Universal Connecting Cable</b> Required with 89955.0000 or 89955.0003.....8564.851
<b>PowerDrive ART1</b>	
	Universal motor system, set incl. power cable, Can Bus connection cable Technical Features: autom. handle and tool recognition, storage function, user-specific parameters, memory function for tools
	<b>Power supply unit 120 V, 50/60 Hz</b> .....2304.0071
	<b>Double-pedal footswitch</b> for PowerDrive ART1 (Series 2304).....2304.901



# COMBIDRIVE EN – High-speed motor system









High-speed motor system	
<b>Burrs with distal protection</b>	
	<b>Round burr, tungsten carbide</b> , burr $\varnothing$ 3.0 mm, WL 350 mm, single-use, sterile, pack of 5.....40990.3730
	<b>Support sleeve</b> , with distal guard, OD 4.0 mm.....82970.1330
	<b>Round burr, diamond</b> , burr $\varnothing$ 3.0 mm, WL 350 mm, single-use, sterile, pack of 5.....40990.3930
	<b>Support sleeve</b> , with distal guard, OD 4,0 mm.....82970.1330
<b>Burrs without distal protection</b>	
	<b>Round burr, diamond</b> , burr $\varnothing$ 3.7 mm, WL 350 mm, single-use, sterile, pack of 5.....40990.3940
	<b>Support sleeve, open, extra visibility</b> , OD 4.0 mm.....88100.1340
<b>High-speed handpiece</b>	
	<b>Handpiece, angled, with adapter</b> 40,000 rpm, INTRA-interface.....82950.1301
<b>Electronic Motor</b>	
	<b>Electronic motor</b> , medium, with connecting cable.....80951.0002
<b>High-speed motor system COMBIDRIVE EN</b> For use with high-speed accessories and accessories for Power Stick M5 (see page 14)	
	<b>COMBIDRIVE EN Set</b> High-speed motor system incl. control unit, power cable, footswitch, electronic motor, connection cable, angled hand-piece and cleaning accessories.....20951.011




# VERTEBRIS stenosis








## Surgimax – Radio frequency generator




RF unit – Surgitron	
	<b>Bipolar accessories – ablation electrodes</b>
	<b>Bipolar hollow sphere electrode</b> distal head part D 2.9 mm, WL 330 mm, disposable .....899364300
	<b>Bipolar hollow sphere electrode</b> distal head part D 3.4 mm, WL 330 mm, disposable .....899364400
	<b>Electrode handle, bipolar</b> for mounting bipolar hollow sphere electrodes, fixed cable with US 2 PIN connector, cable length 3 m, reusable .....899364100
Bipolar accessories – Trigger-Flex	
	<b>Trigger-Flex bipolar RF probe</b> Trigger-Flex RF probe, 31 cm WL.....40990.0031
Surgi-Max® Plus RF Generator	
	<b>Radio frequency generator</b> Surgi-Max® Plus Dual 4 MHz radio frequency unit, power monopolar max. 120 W / 4 MHz, bipolar max. 120 Watt / 1.7 MHz, incl. foot-switch (2345.901) and power cable (7101.54).....2345.0011

## Consumables and replacement parts

for TipControl®:	
	<b>Burr insert, sterile</b> (pack of 5) ..... 499751704
	<b>Wrench</b> ..... 15372005
	<b>Irrigation adapter, complete (M5)</b> ..... 15261106

	<b>Sealing cap attachment</b> incl. 10 sealing caps (89.00) ..... 8792.452
	<b>Sealing caps</b> opening 0.75 mm for instruments up to D 2.4 mm, black, pack of 10..... 89.00
	<b>Sealing caps</b> opening 2.7 mm for instruments over 3.4 to 5.1 mm, blue, pack of 10..... 89.02
	<b>Sealing membrane</b> ..... 15479.006
	<b>Irrigation adapter</b> complete, rotatable..... 15461.034
	<b>O-rings</b> for irrigation toggle (15461.034) pack of 10 ..... 9500.113
	<b>Fog reduction agent</b> for endoscopes, disposable, pack of 10 ..... 102.02

### Cleaning brushes

<b>Cleaning brushes, single-use, pack of 10</b>																																	
	<table border="1"> <thead> <tr> <th>Brush Diameter (mm)</th> <th>Brush Length (mm)</th> <th>Overall Length (mm)</th> <th>Item</th> </tr> </thead> <tbody> <tr> <td>11</td> <td>100</td> <td>600</td> <td>40601</td> </tr> <tr> <td>5</td> <td>75</td> <td>600</td> <td>40605</td> </tr> <tr> <td>5</td> <td>50</td> <td>305</td> <td>40606</td> </tr> <tr> <td>5</td> <td>10</td> <td>250</td> <td>468691</td> </tr> <tr> <td>2</td> <td>15</td> <td>1200</td> <td>7990001</td> </tr> <tr> <td>3</td> <td>20</td> <td>1200</td> <td>7990002</td> </tr> <tr> <td>.85</td> <td>10</td> <td>1200</td> <td>7990003</td> </tr> </tbody> </table>	Brush Diameter (mm)	Brush Length (mm)	Overall Length (mm)	Item	11	100	600	40601	5	75	600	40605	5	50	305	40606	5	10	250	468691	2	15	1200	7990001	3	20	1200	7990002	.85	10	1200	7990003
	Brush Diameter (mm)	Brush Length (mm)	Overall Length (mm)	Item																													
	11	100	600	40601																													
	5	75	600	40605																													
	5	50	305	40606																													
	5	10	250	468691																													
	2	15	1200	7990001																													
3	20	1200	7990002																														
.85	10	1200	7990003																														
Utility brush for external surfaces, disposable, 3/pkg	40999012																																
<b>O-rings</b> for irrigation adapter 892209310 ..... 15364285																																	

# Literature

**Ruetten S, Komp M, Hahn P, Oezdemir S.**

Decompression of lumbar lateral spinal stenosis: full-endoscopic, interlaminar technique. *Oper Orthop Traumatol* 2013;DOI 10.1007/s00064-012-0195-2

**Ruetten S, Komp M, Hahn P, Oezdemir S.**

In: Haertl R, Korge A (eds) *Minimally Invasive Spine Surgery – Techniques, Evidence, and Controversies*. Thieme, Stuttgart 2012, pp 59-62

**Ruetten S.**

Equipment for full-endoscopic spinal surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 59-62

**Ruetten S.**

Endoscopic lumbar disc surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 303-308

**Ruetten S.**

Full-endoscopic operations of the spine in disk herniations and spinal stenosis. *Surg Technol Int* 2011;XXI:284-298

**Ruetten S.**

Full-endoscopic interlaminar lumbar discectomy and spinal decompression. In: Kim DH, Kim K-H, Kim Y-C (eds) *Minimally Invasive Percutaneous Spinal Techniques*. Elsevier, Philadelphia, 2011, pp 351-359

**Komp M, Hahn P, Merk H, Godolias G, Ruetten S.**

Bilateral operation of lumbar degenerative central spinal stenosis in full-endoscopic interlaminar technique with unilateral approach: Prospective 2-year results of 74 patients. *J Spinal Disord Tech* 2011; 24:281-287

**Ruetten S, Komp M, Merk H, Godolias G.**

Surgical treatment for lumbar lateral recess stenosis with the full-endoscopic interlaminar approach versus conventional microsurgical technique: A prospective, randomized, controlled study. *J Neurosurg Spine* 2009;10:476-485

**Ruetten S, Komp M, Merk H, Godolias G.**

Recurrent lumbar disc herniation following conventional discectomy: A prospective, randomized study comparing full-endoscopic interlaminar and transforaminal versus microsurgical revision. *J Spinal Disord Tech* 2009;22:122-129

**Ruetten S.**

Vollendoskopische Operationen der Lendenwirbelsäule. In: Jerosch J, Steinleitner W (ed) *Minimal invasive Wirbelsäulen-Intervention*. Deutscher Ärzte Verlag, Köln, 2009, pp 515-528

**Ruetten S, Komp M, Merk H, Godolias G.**

Full-endoscopic anterior decompression versus conventional anterior decompression and fusion in cervical disc herniations. *Int Orthop* 2008;33:1677, DOI 10.1007/s00264-008-0684-y

**Ruetten S, Komp M, Merk H, Godolias G.**

Full-endoscopic cervical posterior foraminotomy for the operation of lateral disc herniations using 5.9-mm endoscopes: A prospective, randomized, controlled study. *Spine* 2008;33:940-948

**Ruetten S, Komp M, Merk H, Godolias G.**

Full-endoscopic interlaminar and transforaminal lumbar discectomy versus conventional microsurgical technique: A prospective, randomized, controlled study. *Spine* 2008;33:931-939

**Ruetten S, Komp M, Merk H, Godolias G.**

A new full-endoscopic technique for cervical posterior foraminotomy in the treatment of lateral disc herniations using 6.9-mm endoscopes: prospective 2-year results of 87 patients. *Minim Invas Neurosur* 2007; 50:219-226

**Ruetten S, Komp M, Merk H, Godolias G.**

Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disc herniations via the interlaminar and lateral transforaminal approach. *J Neurosurg Spine* 2007;6:521-530

**Ruetten S, Komp M, Godolias G.**

A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6 mm endoscopes: Prospective 2-year results of 331 patients. *Minim Invasive Neurosur* 2006;49:80-87

**Ruetten S, Komp M, Godolias G.**

An extreme lateral access for the surgery of lumbar disc herniations inside the spinal canal using the full-endoscopic uniportal transforaminal approach. – Technique and prospective results of 463 patients. *Spine* 2005;30:2570-2578

**Ruetten S.**

The full-endoscopic interlaminar approach for lumbar disc herniations. In: Mayer HM (ed) *Minimally Invasive Spine Surgery*. Springer, Berlin Heidelberg New York, 2005, pp 346-355





Markkinoija Suomessa



**Mediq Suomi Oy**  
Riihitontuntie 7 D, 02200 Espoo  
puh. 020 112 1510  
asiakaspalvelu@mediq.com  
www.mediq.fi