# **RAPID DIAGNOSTICS**

# OSON<sup>®</sup> Strep A Test

FOR THE DETECTION OF STREPTOCOCCUS A. ANTIGEN

# FOR USE WITH THROAT SWABS



# FEATURES:

- ✓ VERY ACCURATE; 96% sensitivity; 98% specificity versus culture
- ✓ FAST; Results in 5 minutes
- ✓ QC INSIDE; Includes 2 additional test sticks for external QC testing
- ✓ EASY-TO-READ; Two-color results
- WORLD-CLASS TECHNICAL SUPPORT; Staffed with certified medical technologists and skilled professionals

141 OR 141E (CE MARKED)

141E-20 (CE MARKED)

### Place Read results at 10x test stick **5 MINUTES** into the Insert swab and mix. extracted POSITIVE sample. ぼ( )) Strep A Strep A Strep A Add 3 drops Let stand NEGATIVE Reagent 1 and 1 MINUTE. Remove swab 3 drops Reagent 2. 2 Strep A Strep A Strep A from test tube. **ORDERING INFORMATION** CONFIGURATION CATALOG NUMBER

# OSOM<sup>®</sup> STREP A TEST OSOM<sup>®</sup> STREP A TEST

TEST PROCEDURE / RESULTS



Mediq Suomi Oy Luomanportti 3, PL 115 02201 Espoo, puh. 020 112 1500 www.mediq.fi

# **SEKISUI** DIAGNOSTICS

## INTERNATIONAL

**50 TESTS** 

20 TESTS

Sekisui Diagnostics (UK) Limited 50 Gibson Drive Kings Hill, West Malling, Kent, ME19 4AF, UK Phone: +44 (0) 1732 220022 Fax: +44 (0) 1732 220024 Email: info@sekisuidiagnostics.com ISO 9001:2000 Certified, ISO 13485:2003 Certified



OSOM®, the OSOM Training logo and the QC Inside logo are registered trademarks of Sekisui Diagnostics, LLC. ©2011 Sekisui Diagnostics, LLC – All rights reserved.