

M Biliary Stent



Smooth Side-by-Side Stenting Procedure
with A New 6Fr Delivery System

HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

Niti-S™ M Biliary Stent

for malignant hilar strictures



Features

- Braided and a weaving construction
 - Foreshortening reduced while preserving optimal conformability
- Smooth Side-by-Side Stenting Procedure
 - 2 (Two) 6Fr delivery systems can be introduced simultaneously into the working channel for the side-by-side stenting procedure at the hilar biliary stricture
- Radiopaque markers: 3 (three) at both ends & 2 (two) in the middle

Ordering Information

Endoscopic Approach							
Code	Stent		Delivery				
	Diameter (mm) **	Length (cm) ##	Profile (fr)	Usable Length (cm)			
BN0604-6	6	4	6	180			
BN0605-6		5					
BN0606-6		6					
BN0607-6		7					
BN0608-6		8					
BN0609-6		9					
BN0610-6		10					
BN0612-6		12					
BN0804-6		8			4	6	180
BN0805-6					5		
BN0806-6	6						
BN0807-6	7						
BN0808-6	8						
BN0809-6	9						
BN0810-6	10						
BN0812-6	12						

Percutaneous Approach							
Code	Stent		Delivery				
	Diameter (mm) **	Length (cm) ##	Profile (fr)	Usable Length (cm)			
TN0604-6	6	4	6	50			
TN0605-6		5					
TN0606-6		6					
TN0607-6		7					
TN0608-6		8					
TN0609-6		9					
TN0610-6		10					
TN0612-6		12					
TN0804-6		8			4	6	50
TN0805-6					5		
TN0806-6	6						
TN0807-6	7						
TN0808-6	8						
TN0809-6	9						
TN0810-6	10						
TN0812-6	12						

*7Fr delivery systems are also available