

andorate

Precision, Reliability, Consistency



GA Health Story

In 1995 in the midst of the Variant Creutzfeldt-Jakob disease (vCJD) outbreak in the UK, Ken Mc Cabe, the co-founder of GA Health, was invited by Manchester Royal Eye Hospital to develop a solution to prevent the potential transmission of this terrible disease during eye surgery. He took on the challenge and succeeded in being the pioneer in developing ophthalmic instruments for single-patient use in these critical procedures, forever changing infection prevention in ophthalmology.

Today, GA Health is a trusted name for infection prevention especially in endoscopy.

At GA Health, we don't just design products; we solve problems.

It's why medical professionals worldwide count on us, day in, day out.

We bring years of experience, innovative thinking, and a fresh approach to each Infection prevention product we develop.

This is matched only by our passion for the best patient outcomes.

Each solution is thoughtfully designed using innovative technology. And we consistently strive to provide greener, more sustainable solutions.

With our global operation, strong supply chain, and international presence, we continue to be able to offer the reliability and efficiency that our customers need.

Wherever we are, we will never lose our core values that drive us towards innovation and better patient outcomes.

Values

				
Passion	Quality	Integrity	Customer Focus	Simplicity



GA Health is a company with a mission: to develop effective solutions for infection prevention for patients through innovative design that prioritises our users' needs.

We are Global-Asia Health, spanning the globe and providing infection prevention solutions to hospitals and daycare centres in United States, United Kingdom, Europe and Australasia.

We bring innovative design and a dedication to produce consistent and reliable products that benefit both patients and health professionals. Together, we can make GA Health the world's leading solution provider for infection prevention.



Early screening can significantly improve the gastrointestinal disease survival rate. To meet the increasing needs for infection prevention in gastroenterology in Europe, United Kingdom, United States and Australasia, GA Health proudly presents Andorate, a portfolio of high quality and reliable devices for infection prevention from endoscopy preparation, procedure to endoscope cleaning, transport and storage.

Andorate's vision is to serve in every step of the endoscopy cycle.



Ken McCabe
Co-founder & President



Cindy Ye
Co-founder & CEO

Sustainability

In early 2020 GA Health began looking at ways to reduce the packaging we use. We have successfully completed the first phase of this project.

We have modified the layout of our Instructions For Use (IFU). IFUs contain essential information on directions for use of devices. The new layout is more compact reducing our paper usage by 366%.

The introduction of new packaging configurations has also reduced shipping volumes by up to 37%.

Vegetable ink will soon be applied in printing on the packaging.

In addition to these improvements, we use the paper which is Forest Stewardship Council® (FSC®) certified. This gives us the reassurance that the paper we use is sustainably sourced.

These new measures are just the beginning. We continue to research and implement changes focusing on sustainability in our operations and products.

Contents

- 05 Andorate on Olympus® video scope system
- 06 Andorate on Fujifilm® video scope system
- 07 Andorate on Pentax® video scope system

Endo Procedure

- 10 Valves Set
- 13 Biopsy Valve
- 15 Air/Water Bottle Tubing
- 18 Irrigation Connector
- 20 Irrigation Tubing
- 21 Air Supply
- 23 Polyp Trap
- 25 Bite Block

Endo Cleaning

- 28 Cleaning Adapter

Endo Transport & Storage

- 31 Tip Guard

Ordering Information

- 33 All products

Andorate on Olympus® video scope system



Gentle and effective suction performance

Suction Valve
GAR004A



Sufficient insufflation;
Effective lens cleaning

Air/ Water Valve
GAR004B



Soft seal for maintaining the pressure and provide the passage for instrument

Biopsy Valve
GAR028-S




One set for 24-hour air/water delivery

Air/Water Bottle Tubing with CO2 Extend
GAR025CO2EXT



Reusable stainless steel adapter

CO2 Source Tubing Adapter
GAR093



Easy lock and easy open chambers

ThomasTrap™ Polyp Trap with Two Removable Chambers
GAR040



Prevents back-flow during irrigation

Auxiliary Water Connector
GAR027



Stable and reliable irrigation

Irrigation Tubing
GAR083

Andorate on Fujifilm® video scope system



Gentle and effective suction performance

Suction Valve
GAR087A



Sufficient insufflation;
Effective lens cleaning

Air/ Water Valve
GAR087B



Soft seal for maintaining the pressure and provide the passage for instrument

Biopsy Valve
GAR028-S



Ideal for quick irrigation through the biopsy channel

Biopsy Irrigator
GAR015




One for 24-hour air/water delivery

Air/Water Bottle Tubing with CO2 Extend
GAR090CO2EXT



Reusable stainless steel connector

CO2 Source Tubing Adapter
GAR108



Click to open and close

PeralCatch™ Polyp Trap with one fixed chamber
GAR014



Prevents back-flow during irrigation

Auxiliary Water Connector
GAR064



Stable and reliable irrigation

Irrigation Tubing
GAR083

Andorate on Pentax® video scope system



Gentle and effective
suction performance

Suction Valve
GAR037A



Sufficient insufflation;
Effective lens cleaning

Air/ Water Valve
GAR037B



Soft seal for maintaining
the pressure and provide
the passage for instrument

Biopsy Valve
GAR018-S




Prevents back-flow
during irrigation

Auxiliary Water Connector
GAR048



One set for 24-hour
air/water delivery

Air/Water Bottle Tubing
GAR076



Double filters,
click to open and release

Click4Catch™ Polyp Trap
with Four Removable Chambers
GAR082



Stable and reliable
irrigation

Irrigation Tubing
GAR083

Endo Procedure

Endo Procedure Valve Set

Why Andorate valves perform so well?

The products are precisely designed and made with as tight as **0.005mm** tolerance to match the endoscope channels. 0.005mm is only 10% of the average human hair. To achieve such high precision, the hard work is on manufacturability design, material selection, temperature and shrinkage control.

People are the most important reason why GA Health can make such high-quality products. We start from observing and getting to understand how the doctors and nurses complete an endoscopy from patient preparation, procedure, patient discharge to scope reprocessing.

Our passion in designing solutions that can improve the productivity and safety on the field drives all of us along. We only produce the products that we are happy to use on our loved ones.

Three basic notes to identify a good suction and air water valve

1. Action

When the valves are in use, they should be able to quickly and stably rebound upon finger releases. Weak or uneven rebound can cause stickiness, clogging and even jamming during the procedure.

2. Appearance

The edges should be smooth without any jaggedness. Rough edge and surface can cause stickiness, clogging and even jamming during the procedure.

3. Strength

The valves should be able to mount on the port steadily. During the procedure, there is a lot pressure inside the channels. Valves with weak attachment can cause disruption to the procedure in different forms.

Coming in various configurations, Andorate valves are compatible with Olympus, Pentax, and Fujifilm GI endoscopes.

Valve Set

Compatible with Olympus® GI endoscope



Basic Valve Set

Suction, air/water, biopsy valves
3-piece set
GAR004

For Olympus®
140/160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 50 units



Basic Valve Set with Irrigation Connector

Suction, air/water, biopsy valves, irrigation connector
4-piece set
GAR046

For Olympus®
160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 50 units



Basic Valve Set with Air/water Cleaning Adapter

Suction, air/water, biopsy valves, air/water cleaning adapter
4-piece set
GAR3804

For Olympus®
140/160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 10 units



Basic Valve Set with Irrigation Connector, Air/water Cleaning Adapter

Suction, air/water, biopsy valves, irrigation connector, air/water cleaning adapter
5-piece set
GAR3846

For Olympus®
160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 10 units



Combo Valve Set

Suction, air/water, biopsy valves, irrigation connector, air/water cleaning adapter, suction channel cleaning adapter
6-piece set
GAR387546

For Olympus®
160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 10 units

Compatible with
Pentax® GI endoscope



Basic Valve Set

Suction, air/water, biopsy valves
3-piece set
GAR037

For Pentax®
90/i10/i10c Series GI Endoscopes

Qty. per box 50 units



Basic Valve Set
with Irrigation Connector

Suction, air/water, biopsy valves,
irrigation connector
4-piece set
GAR081

For Pentax®
90/i10/i10c Series GI Endoscopes

Qty. per box 50 units

Compatible with
Fujifilm® GI endoscope



Basic Valve Set

Suction, air/water, biopsy valves
3-piece set
GAR087

For Fujifilm®
700 Series GI Endoscopes

Qty. per box 50 units



Basic Valve Set
with Irrigation Connector

Suction, air/water and biopsy valves
and auxiliary water connector
4-piece set
GAR092

For Fujifilm®
700 Series GI Endoscopes

Qty. per box 50 units



Basic Valve Set
with Irrigation Connector,
Air/water Cleaning Adapter

Suction, air/water and biopsy valves
and auxiliary water connector
5-piece set
GAR105092

For Fujifilm®
700 Series GI Endoscopes

Qty. per box 10 units

Endo Procedure Biopsy Valve

The biopsy valve is the most underrated part around the scope. FDA classifies biopsy valve as the Class 2 product which gives the biopsy valve the attention it deserves.

Why the biopsy valve is so important?

- It maintains the air pressure during the procedure.
- It seals off the biohazards inside the endoscope from leaking.
- It provides the passage for the biopsy instruments to get through the endoscope into the patient.

Three basic notes to identify a workable biopsy valve

1. Feel

The biopsy valve should feel not too soft or too hard. Both circumstances can lead to the loss of air inside the working channels. Hard biopsy valve poses a difficulty in inserting instrument, which could cause damage to the instrument.

2. Appearance

The edges should be smooth without any jaggedness. The biopsy valve has indirect contact with the patient by offering a passage for the instrument to get inside the patient. Jaggedness or rough finishes may be introduced into the channel or even to the patient.

3. Strength

The biopsy valve places a critical role as a bio-bacterial seal. It is completely unacceptable that the biopsy valve is not mounted firmly enough on the port. When the biopsy valve pops off the port during a procedure, the contaminated materials from the working channels will immediately get out of the scope. This causes contamination to the patient, to the health professionals and to the procedure room.

The lid of the biopsy valve also presents another risk of cross contamination if there isn't enough strength. It can pop open or leak during the procedure.

Using the right biopsy valve a basic step for an effective and safe endoscopy.

Biopsy Valve



Biopsy Valve

GAR018

For Pentax® GI Endoscopes

NON-STERILE

Qty. per box 100 units

GAR018-S

For Pentax® GI Endoscopes

STERILE

Qty. per box 100 units



Biopsy Valve

GAR028

For Fujifilm® & Olympus® GI Endoscopes

NON-STERILE

Qty. per box 100 units

GAR028-S

For Fujifilm® & Olympus® GI Endoscopes

STERILE

Qty. per box 100 units

Endo Procedure Air/Water Bottle Tubing

The endoscope can't create the air and water system without the air/ water bottle tubing. The air water bottle tubing works with the air/water valve to insufflate the intestine and clean the lens for better view.

Three basic notes to identify a workable air/water bottle tubing

1. The softness

The tubing should be soft enough, to accommodate between the water bottle and the scope during the 24-hour use.

2. The strength


For example, there is up to **80 psi** pressure during a colonoscopy. The connector must be able to sustain the pressure and stay on the scope all the time during the endoscopy. Weak connection point could cause disruption to the procedure

3. The compatibility

The water bottle thread may vary from market to market, brand to brand. Intelligent engineering is needed to make a cap that is compatible with different water bottles. If there is insufficient insufflation, the compatibility of the cap to the water bottle should be on the check-list.

Air/Water Bottle Tubing

For Olympus® GI endoscope



Air/Water Bottle Tubing
GAR025

For Olympus®
140/160/180/190/240/260/290 Series
GI Endoscopes

Qty. per box 20 units



Air/Water Bottle Tubing
with CO₂ Extend
GAR025CO2EXT

For Olympus®
140/160/180/190/240/260/290 Series
GI Endoscopes

Qty. per box 10 units

For Pentax® GI endoscope



Air/Water Bottle Tubing
GAR072

For Pentax® 90 Series
GI Endoscopes

Qty. per box 20 units



Air/Water Bottle Tubing
with CO₂ Extend
GAR072CO2EXT

For Pentax® 90 Series
GI Endoscopes

Qty. per box 10 units

Air/Water Bottle Tubing

For Fujifilm® GI endoscope



Air/Water Bottle Tubing
GAR088

For Fujifilm® 500/600 Series
GI Endoscopes

Qty. per box 20 units



Air/Water Bottle Tubing
with CO₂ Extend
GAR088CO2EXT

For Fujifilm® 500/600 Series
GI Endoscopes

Qty. per box 10 units



Air/Water Bottle Tubing
GAR090

For Fujifilm® 700 Series
GI Endoscopes

Qty. per box 20 units



Air/Water Bottle Tubing
with CO₂ Extend
GAR090CO2EXT

For Fujifilm® 700 Series
GI Endoscopes

Qty. per box 10 units

Endo Procedure Irrigation Connector

What is the most important job of an irrigation connector besides connecting the irrigation tubing to the scope?

That is to prevent the contaminated water in the endoscope channel from backflowing to the 24-hour multi-patient use water bottle.

Too much resistance will reduce the irrigation water flow. Too little resistance will fail to prevent the back flow. Andorate knows the key to this balance.

Irrigation Tubing

What can possibly go wrong with the irrigation during the GI procedure?

Dripping, blowing up, breaking-up, back-flow, no priming could appear.

Andorate irrigation tubing gives a 24-hour stable and reliable flow. With the anti back-flow valve, the 24-hour use irrigation tubing with the single-patient use connector provides a safe GI irrigation.

Use Andorate single-patient use connector with the 24-hour use irrigation tubing to provide a safe GI irrigation.

Irrigation Connector



Biopsy Irrigator

GAR015

For Fujifilm® & Olympus®
GI Endoscopes

Qty. per box 100 units



Auxiliary Water Connector

GAR027

For Olympus®
160/180/190/240/260/290
Series GI Endoscopes

Qty. per box 100 units



Auxiliary Water Connector

GAR048

For Pentax® 90/i10/i10c Series
GI Endoscopes

Qty. per box 100 units



Auxiliary Water Connector

GAR064

For Fujifilm® 500/600/700 Series
GI Endoscopes

Qty. per box 100 units

Irrigation Tubing



Irrigation Tubing

GAR073

For Olympus® OFP and
Endo Stratus™ EGA-500 Irrigation Pump

Qty. per box 10 units



Irrigation Tubing

GAR083

For Olympus® OFP, Endo Stratus™
EGA-500 and EGP-100 Irrigation Pump

Qty. per box 10 units

Endo Procedure Air Supply

Carbon dioxide (CO₂) is widely used for insufflation in endoscopy.

CO₂ insufflation is demonstrated in the literature to cause less abdominal discomfort as it is quickly reabsorbed by the body. Andorate presents a very high quality of CO₂ source tubing with very fine micro filter, which particularly works for the insufflator with a pressure sensing mode.

Air Supply



CO₂ Source
Tubing Adapter

GAR093

For Olympus®
UCR CO₂ Insufflator

Qty. per box 100 units



CO₂ Source
Tubing Adapter

GAR108

For Fujifilm®
GW-100 Insufflator

Qty. per box 100 units



CO₂ Source Tubing with
Luer Connector

GAR047

Qty. per box 10 units

Endo Procedure Poly Trap

How to make better use of the traditional polyp trap in the GI procedure?

What if the nurse can see inside the trap at every angle? What if the trap could give clear indication to the nurse in opening, twisting, closing or separating the samples? What if there is extra filter to ensure the sample has no way to escape?

Andorate polyp trap range is designed with all these benefits in mind. From single chamber to four chambers, there must be one for your need. If not, contact us any time. We would love to hear your feedback.

Polyp Trap



PearlCatch™ Polyp Trap
with one fixed chamber
GAR014

Qty. per box 50 units



ThomasTrap™ Polyp Trap
with two removable chambers
GAR040

Qty. per box 20 units



Click4Catch™ Polyp Trap
with four removable chambers
GAR082

Qty. per box 12 units



Basic Polyp Trap
with four fixed chambers
GAR023

Qty. per box 12 units

Endo Procedure Bite Block

What can possibly go wrong with a bite block?

On FDA's Maude Adverse Event Report, there are dozens of reports related to bite blocks such as laceration on patient's tongue, lip and face, or a crack happens when it is bitten too hard.

To ensure structure integrity achieve smooth edges, we put extra effort on developing its manufacturability. Our bite blocks are ergonomically designed and provide extra comfort to the patients with soft comfortable silicone rubber

Bite Block



Bite Block
66Fr with strap, extra space for finger
SDPC0210

Qty. per box 100 units



Bite Block
66Fr with strap and tooth protector,
extra space for finger
SDPC0210F

Qty. per box 100 units



Bite Block
55.5Fr with strap
SDPC0209A

Qty. per box 100 units



Bite Block
55.5Fr with strap and tooth protector
SDPC0209

Qty. per box 100 units



Bite Block
66Fr with oxygen connection
GAR044

Qty. per box 100 units



Bite Block
66Fr with oxygen connection
and tooth protector
GAR044F

Qty. per box 100 units

Endo Cleaning

Endo Cleaning Cleaning Adapter

To improve the productivity in endoscopy is the core objective of Andorate.

The cleaning adapters, designed and manufactured at a high engineering level, are our best sellers.

They come in different configurations to offer convenience at the point-of-use cleaning.

Cleaning Adapter



For Olympus®
GI Endoscope



For Fujifilm®
GI Endoscope



Air/ Water Channel
Cleaning Adapter
GAR038

For Olympus® GI Endoscopes

Qty. per box 100 units



Air/Water Channel
Cleaning Adapter
GAR105

For Fujifilm® 700 Series GI Endoscopes

Qty. per box 100 units



Basic Cleaning
Adapter Set

Suction & air/water channel cleaning adapter
GAR3875-NS

For Olympus® GI Endoscopes

Qty. per box 50 units

Endo Transport & Storage

EndoTransport & Storage Tip Guard

How to reduce the wear-and-tear during the endoscope transport and storage?

Ordinary products like the sponge do not allow the scope to dry completely, and the mesh tube does not provide strong protection. When we were putting the ideas together, we set a clear objective of this product is to protect the scope during the transport and storage with zero moisture.


Zero moisture requires great ventilation structure.

The greater ventilation, the less contact with the scope.

Andorate tip guard is designed to softly attached to the scope at all time during the transport and storage with as little as possible contact.

We even gave ourselves a further push to develop the tip guard to fit most scopes.

Tip Guard



Secure Attachment without Moisture

For GI Endoscope

Effective Moisture Ventilation

Minimum contact with maximised ventilation to prevent bioburden growing.

Secure Attachment

Silicon discs allow secure attachment to the scope avoiding possible breakage during transport and storage

One Size Solution

One size fits most scopes to provide easy supply management for the users.

Robust Shock Absorption

Designed to absorb shock and vibration effectively with the rib structure.

Endoscope Tip Guard

GAR107
NON-STERILE

Qty. per box 50 units

GAR107-S
STERILE






Qty. per box 50 units

Ordering Information

The Valve Set Configurations

for Olympus® GI Endoscope

Suction Valve Air/ Water Valve Biopsy Valve Auxiliary Water Connector Air/Water Channel Cleaning Adapter Suction Channel Cleaning Adapter

Combination						
GAR004	●	●	●			
GAR046	●	●	●	●		
GAR3804	●	●	●		●	
GAR3846	●	●	●	●	●	
GAR3875-NS					●	●
GAR387504	●	●	●		●	●
GAR387546	●	●	●	●	●	●

The Valve Set Configurations

for Fujifilm® GI Endoscope

Suction Valve Air/ Water Valve Biopsy Valve Auxiliary Water Connector Air/Water Channel Cleaning Adapter

Combination					
GAR087	●	●	●		
GAR092	●	●	●	●	
GAR105087	●	●	●		●
GAR105092	●	●	●	●	●

Endo Procedure

Valves Set

Product Code	Description	Qty Per Box
GAR004C	2-piece Valve Kit (Suction, Air/Water Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR004	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR046	4-piece Valve Kit (Suction, Air/Water, Biopsy Valves & Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR3804	A/W Channel Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR3846	A/W Channel Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water, Biopsy Valves & Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR387504	2-piece Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR387546	2-piece Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water, Biopsy Valves and Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR065	3-piece Valve Kit (Suction-Metal Shaft, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR066	4-piece Valve Kit (Suction-Metal Shaft, Air/Water, Biopsy Valves & Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR037C	2-piece Valve Kit (Suction, Air/Water Valves) For Pentax® 90/i10/i10c Series GI Endoscopes, Sterile	50 units
GAR087C	2-piece Valve Kit (Suction, Air/Water Valves) For Fujifilm® 700 Series GI Endoscopes, Sterile	50 units
GAR087	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Fujifilm® 700 Series GI Endoscopes, Sterile	50 units
GAR092	4-piece Valve Kit (Suction, Air/Water & Biopsy Valves & Auxiliary Water Connector) For Fujifilm® 700 Series GI Endoscopes, Sterile	50 units
GAR105087	A/W Channel Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Fujifilm® 700 Series GI Endoscopes, Sterile	10 units
GAR105092	A/W Channel Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water, Biopsy Valves & Auxiliary Water Connector) For Fujifilm® 700 Series GI Endoscopes, Sterile	10 units
GAR037C	2-piece Valve Kit (Suction, Air/Water Valves) For Pentax® 90/i10/i10c Series GI Endoscopes, Sterile	50 units

Endo Procedure

Valves Set

Product Code	Description	Qty Per Box
GAR037	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Pentax® 90/i10/i10c Series GI Endoscopes, Sterile	50 units
GAR081	4-piece Valve Kit (Suction, Air/Water & Biopsy Valves & Auxiliary Water Connector) For Pentax® 90/i10/i10c Series GI Endoscopes, Sterile	50 units

Biopsy Valve

GAR028	Biopsy Valve For Olympus® and Fujifilm® GI Endoscopes, Non-Sterile	100 units
GAR028-S	Biopsy Valve For Olympus® and Fujifilm® GI Endoscopes, Sterile	100 units
GAR018	Biopsy Valve For Pentax® GI Endoscopes, Non-Sterile	100 units
GAR018-S	Biopsy Valve For Pentax® GI Endoscopes, Sterile	100 units

Air/Water Bottle Tubing

GAR025	Vacutore® Air/Water Bottle Tubing For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	20 units
GAR025CO2EXT	Vacutore® Air/Water Bottle Tubing with CO₂ Extend For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR088	Vacutore® Air/Water Bottle Tubing For Fujifilm® 500/600 Series GI Endoscopes, Sterile	20 units
GAR088CO2EXT	Vacutore® Air/Water Bottle Tubing with CO₂ Extend For Fujifilm® 500/600 Series GI Endoscopes, Sterile	10 units
GAR090	Vacutore® Air/Water Bottle Tubing For Fujifilm® 700 Series GI Endoscopes, Sterile	20 units
GAR090CO2EXT	Vacutore® Air/Water Bottle Tubing with CO₂ Extend For Fujifilm® 700 Series GI Endoscopes, Sterile	10 units
GAR072	Vacutore® Air/Water Bottle Tubing For Pentax® 90 Series GI Endoscopes, Sterile	20 units
GAR072CO2EXT	Vacutore® Air/Water Bottle Tubing with CO₂ Extend For Pentax® 90 Series GI Endoscopes, Sterile	10 units
GAR076	Vacutore® Air/Water Bottle Tubing with Side Port For Pentax® 90 Series GI Endoscopes, Sterile	20 units

Air Supply

Product Code	Description	Qty Per Box
GAR047	Vacutore® CO ₂ Source Tubing with Luer Connector Non-Sterile	10 units
GAR091	Vacutore® CO ₂ Extend Tubing with Luer Connector Non-Sterile	20 units

Reusable Connector

Product Code	Description	Qty Per Box
GAR093	Vacutore® CO ₂ Source Tubing Adapter For Olympus® UCR CO ₂ Insufflator, Non-Sterile	100 units
GAR108	Vacutore® CO ₂ Source Tubing Adapter For Fujifilm® GW-100 Insufflator, Non-Sterile	100 units

Irrigation Connector

Product Code	Description	Qty Per Box
GAR015	Biopsy Irrigator For Olympus® and Fujifilm® GI Endoscopes, Sterile	100 units
GAR027	Aquapulse® Auxiliary Water Connector For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	100 units
GAR064	Aquapulse® Auxiliary Water Connector For Fujifilm® 500/600/700 Series GI Endoscopes, Sterile	100 units
GAR048	Aquapulse® Auxiliary Water Connector For Pentax® 90/i10/i10c Series GI Endoscopes, Sterile	100 units

● Olympus® ● Fujifilm® ● Pentax®

Irrigation Tubing

Product Code	Description	Qty Per Box
GAR073	Aquapulse® Irrigation Tubing For Olympus® OFP and Endo Stratus™ EGA-500 Irrigation Pump, Sterile	10 units
GAR083	Aquapulse® Irrigation Tubing For Olympus® OFP and Endo Stratus™ EGA-500 and Endogator™ EGP-100 Irrigation Pump, Sterile	10 units

Polyp Trap

Product Code	Description	Qty Per Box
GAR014	PearlCatch™ Single Chamber Polyp Trap Non-Sterile	50 units
GAR023	Four Chamber Polyp Trap Non-Sterile	12 units
GAR040	ThomasTrap™ Polyp Trap with two removable Chambers Non-Sterile	20 units
GAR082	Click4Catch™ Click into Place Polyp Trap with four removable Capture Chambers Non-Sterile	12 units

Bite Block

Product Code	Description	Qty Per Box
GAR044	Bite Block 66Fr with Oxygen Connection Non-Sterile	100 units
GAR044F	Bite Block 66Fr with Oxygen Connection and Tooth Protector Non-Sterile	100 units
SDPC0209	Bite Block 55.5Fr with Strap and Tooth Protector Non-Sterile	100 units
SDPC0209A	Bite Block 55.5Fr with Strap Non-Sterile	100 units
SDPC0210	Bite Block 66Fr with Strap Non-Sterile	100 units
SDPC0210F	Bite Block 66Fr with Strap and Tooth Protector Non-Sterile	100 units

Cleaning Adapter

Product Code	Description	Qty Per Box
GAR3875	2-piece Cleaning Adapter Kit (Suction & Air/Water Channel Cleaning Adapter) For Olympus® GI Endoscopes, Sterile	25 units
GAR3875-NS	2-piece Cleaning Adapter Kit (Suction & Air/Water Channel Cleaning Adapter) For Olympus® GI Endoscopes, Non-Sterile	50 units
GAR075	Suction Channel Cleaning Adapter For Olympus® GI Endoscopes, Non-Sterile	50 units
GAR105	Air/Water Channel Cleaning Adapter For Fujifilm® 700 Series GI Endoscopes, Non-Sterile	100 units

● Olympus® ● Fujifilm®



Endo Transport & Storage

Tip Guard

Product Code	Description	Qty Per Box
GAR107	Endoscope Tip Guard Non-Sterile	50 units
GAR107-S	Endoscope Tip Guard Sterile	50 units

andorate

Precision, Reliability, Consistency



Headquarters:

GA Health Company Limited

Unit 18, 21/F, Metropole Square
2 On Yiu Street Shatin, N.T.
Hong Kong

+852 2833 9010
info@GAhealth.com
GAhealth.com

Europe:

GA Health Medical Devices Limited

Forest Park Business Campus
Mullingar, N91 VEOH, Co. Westmeath
Ireland

+353 (0)44 931 1050
inquiry@gahealth.com
GAhealth.com

All trademarks are the property of their respective owners
GAH-SO-GEN-C001 Rev 1.0