BioGuard® air/water & suction valve kits air/water cleaning adapter





NOT for Patient Use ONLY for Pre-Cleaning Air / Water Channel DO NOT REUSE



Endoscopy

STERIS Endoscopy is proud to expand our portfolio of disposable infection prevention joins the **BioGuard® air/water & suction valve kits** designed to support pat

The single-use BioGuard air/water & suction valve kits eliminate the need for manual cleaning and reprocessing.

## The BioGuard air/water & suction valve kits...

- are available in 6 unique kit configurations
- provide consistent performance use a new kit for each procedure
- disposability supports adherence to industry guidelines which require traceability with a specific endoscope<sup>1,2</sup>



<sup>1</sup>ANSI/AAMI ST91: 2015 Flexible and semi-rigid endoscope processing in health care facilities. 20, 23, 25, 31, 32, 33. <sup>2</sup>ESGE-ESGNA Position Statement: 2018 Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGNA).

product number	description	
BX00711780*	BioGuard air/water & suction valves - 2 piece kit	50
BX00711783*	BioGuard air/water & suction valves - 3 piece kit <b>Kit includes:</b> BioShield® biopsy valve (BX00711129)	50
BX00711784*	BioGuard air/water & suction valves - 4 piece kit <b>Kit includes:</b> BioShield® biopsy valve (BX00711129) and Torrent® scope connector (BX00711570)	50
BX00711785	BioGuard valves & air/water cleaning adapter kit	50
BX00711786	BioGuard valves (3pc kit) & air/water cleaning adapter kit	50
BX00711787	BioGuard valves (4pc kit) & air/water cleaning adapter kit	50

\*Each BioGuard air/water & suction valve kit (BX00711780, BX00711783, BX00711784) comes sterile, individually packaged and sealed in tray. Compatible with Olympus GI endoscopes (excluding EUS echoendoscopes that utilize a balloon channel). BioGuard air/water cleaning adapter is non-sterile (BX00711788). on solutions. The new **BioGuard<sup>®</sup> air/water cleaning adapter - Olympus** ient safety and support infection prevention practices.

The new BioGuard air/water cleaning adapter is designed specifically for Olympus GI endoscopes<sup>3</sup>, offering a single-use option that flushes the air/water channel during the bedside pre-cleaning process.

## The BioGuard air/water cleaning adapter...

- is disposable, offering an alternative option to manual reprocessing and the need to be tracked with a specific endoscope
- offers a unique appearance and warning tag to prevent use during patient procedures
- individually packaged to ensure use only during the bedside pre-cleaning process



<sup>3</sup>Compatible with Olympus GI endoscopes, excluding EUS echoendoscopes that utilize a balloon channel.

product nu	umber	description	sterile	units/box
BX00711	1788	BioGuard air/water cleaning adapter - Olympus	no	50

For more information on the BioGuard products, contact STERIS Endoscopy today.





Endoscopy

STERIS 5960 Heisley Road Mentor, OH 44060-1834 – USA Customer Service: 800-548-4873 Fax: 440-639-4495 www.steris.com

Olympus<sup>®</sup> is a registered trademark of Olympus Corporation. ©2019 STERIS Endoscopy or its affiliates. All rights reserved. Unless otherwise indicated, all marks denoted with @ or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS or its affiliates.

