

OSOM[®]

Trichomonas Rapid Test



THE ONLY CLIA-WAIVED RAPID TEST FOR THE DETECTION OF *TRICHOMONAS VAGINALIS*


FOR USE WITH VAGINAL SWABS

FEATURES:



- ✓ **ACCURATE**; 83% sensitivity vs. culture; 95% agreement with culture and wet mount combined
- ✓ **FAST**; Results in 10 minutes or less
- ✓ **OBJECTIVE**; Detects the antigen; removes uncertainty associated with live organism instability
- ✓ **WORLD-CLASS TECHNICAL SUPPORT**; Staffed with certified medical technologists and skilled professionals

TEST PROCEDURE / RESULTS



ORDERING INFORMATION

CONFIGURATION

CATALOG NUMBER

○ OSOM[®] TRICHOMONAS RAPID TEST
TRICHOMONAS POSITIVE CONTROL KIT

25 TESTS

181 OR 181E (CE MARKED)
182



Mediq Suomi Oy
Luomanportti 3, PL 115
02201 Espoo, puh. 020 112 1500
www.mediq.fi



INTERNATIONAL
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling, Kent, ME19 4AF, UK
Phone: +44 (0) 1732 220022
Fax: +44 (0) 1732 220024
Email: info@sekisuidiagnostics.com
ISO 9001:2000 Certified, ISO 13485:2003 Certified



MADE IN THE USA