


INLAY OPTIMA[®]

URETERAL STENT

Maximum Efficiency.
Optimum Results.

FEATURES

- Made of a proprietary polymer blend that results in a biocompatible and ultra smooth stent
- Depending on stent size, stent softens from 29% up to 49% at body temperature to promote greater patient comfort*
- Superior resistance to urine calcium salt accumulation compared to two leading competitive stents*
- Gradual taper at tip of stent promotes smooth insertion and atraumatic navigation around obstructions



*In-Vitro Testing**

INLAY OPTIMA[®] Ureteral Stent

Size	Without Guidewire	Nitinol Guidewire
4.7 Fr/14 cm	788414	786414
4.7 Fr/20 cm	788420	786420
4.7 Fr/24 cm	788424	786424
4.7 Fr/26 cm	788426	786426
4.7 Fr/28 cm	788428	786428
4.7 Fr/30 cm	788430	786430
6 Fr/20 cm	788620	786620
6 Fr/22 cm	788622	786622
6 Fr/24 cm	788624	786624
6 Fr/26 cm	788626	786626
6 Fr/28 cm	788628	786628
6 Fr/30 cm	788630	786630
7 Fr/20 cm	788720	786720
7 Fr/22 cm	788722	786722
7 Fr/24 cm	788724	786724
7 Fr/26 cm	788726	786726
7 Fr/28 cm	788728	786728
7 Fr/30 cm	788730	786730
8 Fr/22 cm	788822	786822
8 Fr/24 cm	788824	786824
8 Fr/26 cm	788826	786826
8 Fr/28 cm	788828	786828
8 Fr/30 cm	788830	786830